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MEDISOLARIS Dental Implant Instruction for Use

1-) DENTAL IMPLANTS

Medisolaris (VENUSCON) Dental Implant and Superstructures produced by MEDISOLARIS Sağlık Hizmetleri San. Tic. Ltd. Ltd. It is and is surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations. VENUSCON Dental Implants are bone level dental implants. The target users for Venüscon Dental Implants are dentists.

Superstructures are the device that used to support crowns, bridges or superstructures in the lower and/or upper jaw, used with Dental Implants in completely edentulous or partially edentulous patients.

It is recommended that Medisolaris products not be used for purposes other than their stated intended use.

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1.e Patient Population : Patients with complete bone development (18 years of age or older), not pregnant, not in the rapid puberty period of craniofacial changes and not during lactation, without hard or soft tissue pathology in the implant area, without active periodontal infection or untreated periodontitis, after periodontal treatment; good oral hygiene, good response to periodontal treatment, non/smoker, systemically healthy patients, no unstable endocrine diseases (such as diabetes), no local inflammation or mucosal diseases such as lichen planus, no osteoporosis (Bone resorption), no immunosuppressive disease, no cancer cases requiring chemotherapy / radiotherapy / anti-resorptive drug use, no acquired or acquired coagulation disorders, no intravenous (IV) based bisphosphonate use / use, no history of alcohol and drug addiction, no history of allergy to titanium and stainless steel, no systemic disease requiring regular medication, no mental development disorders such as mental retardation. Venuscon Implants have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. They have not been tested for heating, displacement or image artefact in the MR environment. Venuscon safety in the MR environment is unknown. Scanning a patient with this device may result in patient injury.

1.f Raw Materials and Ingredients: Implants and superstructures are made of ISO 5832-3 Ti6Al4V ELI / GRADE23, surgical instruments are made of 1.4116, hand instruments are made of 1.4057 Stainless Steel. The surface of Medisolaris Dental Implant is roughened by sandblasting and passivation SBM method consisting of hydroxyapatite and tricalcium phosphate to obtain surface microstructure, osseointegration, Ra values 2.0 - 2.5 µm.

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1.h Contraindications:

<ul style="list-style-type: none"> ✓ Having hard or soft tissue pathology in the İmplant area ✓ Those with local inflammation or mucosal diseases such as lichen planus ✓ People with immunosuppressive disease or on medication ✓ Have a history of radiation therapy to the head or neck area ✓ Used/used intravenous (IV) based bisphosphonate ✓ Patients with recent myocardial infarction and cerebrovascular accident 	<ul style="list-style-type: none"> ✓ Those with active periodontal infection or untreated periodontitis ✓ Patients with systemic problems that cannot be controlled and will negatively affect wound healing ✓ A disease that affects osteoporosis and/or bone metabolism ✓ Having a history of alcohol and drug addiction ✓ Have a history of allergy to titanium ✓ Have any disorder that interferes with the ability to maintain adequate daily oral hygiene or plaque control
---	---

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
1.j Potential Side Effects:


• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis

• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening

1.k Placement Configuration of Implants in the Mouth : Since the tooth itself and its root are small in diameter, small diameter implants are placed in areas where the alveolar crest is thinner. After the implant is placed in the alveolar ridge, there should be at least two mm of bone distance around the implant in bucco-lingual and mesio-distal directions. This two mm distance provides the necessary support to the implant during use. In Venuscon series implants; It is recommended that 3.2 mm diameter implants with a narrow platform be placed in the dark area on the side. It is recommended that the implant size to be used is decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME	PRODUCT CODE	PRODUCT NAME
SLV 3208	Venuscon Implant Ø3.2 mm L:8 mm Platform 3.0mm	SLV 3213	Venuscon Implant Ø3.2 mm L:13 mm Platform 3.0 mm
SLV 3210	Venuscon Implant Ø3.2 mm L:10 mm Platform3.0mm	SLV 3216	Venuscon Implant Ø3.2 mm L:16 mm Platform 3.0 mm
SLV 3211	Venuscon Implant Ø3.2 mm L:11.5 mm Platform3.0 mm		

follows.

Likewise, large diameter implants are used in areas where the tooth and bone diameter is thick, that is, there are teeth with more than one root. A large diameter implant is placed in the area where the lower and upper jaw incisors are located, where there is biting force instead of chewing. In Medisolaris implants; Venüscon Ø3.7/4.3/5.0 implants, which have a regular platform, are recommended to be placed in the dark area on the side. It is recommended that the implant size to be used is decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME	PRODUCT CODE	PRODUCT NAME
SLV 3706	Venuscon Implant Ø3.7mm L:6mm Platform 3.0mm	SLV 4306	Venuscon Implant Ø4.3mm L:6mm Platform 3.0mm
SLV 3708	Venuscon Implant Ø3.7mm L:8mm Platform 3.0mm	SLV 4308	Venuscon Implant Ø4.3mm L:8mm Platform 3.0mm
SLV 3710	Venuscon Implant Ø3.7mmL:11.5mmPlatform3.0mm	SLV 4310	Venuscon Implant Ø4.3mm L:10mm Platform 3.0mm
SLV 3711	Venuscon Implant t Ø3.7mm L:11.5mmPlatform3.0mm	SLV 4311	Venuscon Implant Ø4.3 mm L:11mm Platform 3.0mm
SLV 3713	Venuscon Implant Ø3.7mm L:13mm Platform3.0 mm	SLV 4313	Venuscon Implant Ø4.3mm L:13mm Platform3.0mm
SLV 3716	Venuscon Implant Ø3.7mm L:16mm Platform 3.0mm	SLV 4316	Venuscon Implant Ø4.3mm L:1 mm Platform 3.0mm
SLV 5006	Venuscon Implant Ø5.0mm L:6mm Platform 3.0mm	SLV 5011	Venuscon Implant Ø5.0mm L:11.5mmPlatform3.0mm
SLV 5008	Venuscon Implant Ø5.0mm L:8mm Platform 3.0mm	SLV 5013	Venuscon Implant Ø5.0mm L:13mm Platform3.0mm
SLV 5010	Venuscon Implant Ø5.0mm L:10mm Platform 3.0mm	SLV 5016	Venuscon Implant Ø5.0mm L:16mm Platform3.0mm

1.l Packaging and Sterilization : Medisolaris Dental implants are packaged sterile with Caps. The products are disposable. Gamma radiation is used to sterilize implants after packaging. Thanks to sealed packaging, implants are prevented from being affected by external environments. For proper storage conditions, the temperature should be between 18-28°C and the humidity should be between 40-60%. Sterilization is valid until the expiry time stated on the label on the package. Implants in damaged sterile packages should not be used. Contaminated implants cannot be resterilized. Reuse of implants may result in clinical failure and complications including, but not limited to, infection, failure to osseointegrate implants, bone loss, and damage to adjacent structures. The implant surface should not be touched by hand, as this may deteriorate its surface properties. All other superstructures and hand instruments should be sterilized by autoclave. For recommended autoclave

sterilization parameters, see the Superstructure section and *MKK.02 Surgical Instruction for Use*. Do not use the product whose packaging is damaged.

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For detailed information see *MKK.02 Surgical Instruction for Use*. *Recommended used at 700 RPM.

1.q Healing Process of Medisolaris Dental Implants: In general, implants should be left to heal for 2 to 4 months, depending on the patient's bone quality, type and general health status. Immediate loading can be performed with the right case selection and adherence to the immediate loading protocols accepted in the literature. If immediate loading is to be performed, the occlusion environment should be checked. If there are multiple implants, they should be splinted and in cases of complete edentulism, at least 4 implants in the mandible and 6 implants in the maxilla should be connected to each other with a prosthesis. Depending on the bone anatomy, the appropriate number, diameter and length of implants should be placed.

1.r Information to be communicated to the patient after implant placement

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1. After the dental implant, the dentist may prescribe antibiotics and antiseptic mouthwash to the patient if deemed appropriate.
2. A sterile sponge is placed on the area where dental implant surgery is performed to control bleeding. The sterile sponge should be bitten by the patient with light force for 20 minutes. At the end of 20 minutes, the sponge should be removed.
3. The patient should not eat or drink anything for two hours after the operation. The patient should not spit or rinse his/her mouth with any liquid for two hours. After two hours, the patient should not eat hard and hot foods on the first day. Hot foods increase both bleeding and the risk of infection.
4. It is normal to have bleeding in the form of leakage during the first 24 hours after the procedure. To control bleeding, the patient bites gauze for 30 minutes. If the bleeding continues for more than 2 days, you should consult our physician. Aspirin and aspirin compounds should not be taken as painkillers because they cause bleeding tendency.

5. Swelling may occur after implant surgery. This is a normal condition. If the pain and swelling do not disappear after 5-6 days, you should consult your physician.
6. The patient should avoid tampering with the area where the dental implant is placed with the tongue, finger or hard instruments.
7. The patient should rest a little after the implant surgery and avoid strenuous exercise for 48 hours.
8. If the patient is using removable dentures, he/she should not use them until the date that the dentist deems appropriate.
9. In dental implant surgery, the wound surfaces are closed with sutures. The sutures must be removed by the dentist after one week at the earliest. The patient is informed by the dentist about the early period in which the treated tooth part can be used.
10. After dental implant surgery, the patient should continue oral care after the first 24 hours. The patient should use antiseptic mouthwash prescribed by the dentist three times a day for the first week after the surgery. Antiseptic mouthwash should not be used after one week. After antiseptic mouthwash, the patient can gargle with warm salt water. Add 1 teaspoon of salt to 1 glass of water and gargle 3-4 times a day to disinfect the area. This will help clean the healing area.
11. After 48 hours, the surgical site should be lightly brushed. Other teeth can be brushed and cleaned normally after the operation. There should be no food residue in the surgical site and the patient should rinse the mouth with warm water after meals.
12. Smoking delays the healing process and prevents dental implants from staying in the mouth for a long time, so the patient should stop smoking cigarettes and similar tobacco before starting dental implant treatment.
13. It is very important to follow all these instructions mentioned above for the success of dental implant surgery. If there are signs of infection after implant surgery or if the patient experiences nausea, fever or vomiting, he/she should consult his/her dentist.

1.s Long-term care of dental implants: In order to be able to use implants for a lifetime without any problems, it is important that they are properly cared for as described above.

1.t The importance of oral hygiene in dental implant patients: Implants are more susceptible to plaque and calculus build-up, inflammation and bone loss than natural teeth. In implant patients, good oral hygiene is very important. The patient should brush teeth correctly and clean the interdental and interimplant areas correctly. Oral hygiene kits are available for the private use of implant patients, which can be used by the patient. If the patient follows the oral hygiene instructions, inflammation of the periimplant tissues will be minimal. Removal of plaque and calculus deposits from the teeth will reduce the risk of periimplantitis, the most common cause of implant failure.

1.u Regular dental visits and examinations for implant patients: In patients treated with dental implant prostheses, check-ups are usually recommended every 3 months for the first 1-2 years and then every 6 months thereafter. At each visit, oral hygiene instructions are reinforced and the patient's ability to maintain good oral hygiene is checked. Patients with implant-retained removable prostheses are given some special instructions. The patient is told how to insert and remove the prosthesis without causing any damage to the implant or prosthesis. Excessive pressure should not be applied to the restoration placed on the implant. Restorations on implants are not as strong as natural teeth. The patient should avoid eating hard foods that can break the restoration. Avoid irritation of the implant area.

1.v Lifespan of dental implants: The lifespan of dental implants depends on many factors, such as patient health and the care provided by the patient. The lifespan of implants depends on how old the patient had the implants, how they live and how they care for them. For the success of the implants, the patient should maintain good oral hygiene and visit the dentist regularly for check-ups. Although a clear lifetime cannot be defined for the device, it is thought to have a lifetime of 10 years and longer under normal conditions. Implants have a shelf life of 5 years and superstructures do not have a shelf life.

Implants are an effective solution for missing teeth and can be used for a lifetime with good care. Brushing, flossing and regular visits to the dentist after the implant prosthesis are important to prolong the implant's stay in the mouth. Implants can replace a single tooth, a group of teeth or all upper and lower jaw teeth.

1.w Opening Sterile Packages: To take Medisolaris Dental implants from the product packaging with an angldurva, implant keys suitable for angldurva are used. To take the implants manually, the implant key suitable for manual removal is used. When taking the implants, the size suitable for the area to be placed in the mouth should be selected. For all diameters of the Venuscon implant, the 2.3 mm hexagon crimping spanner included in the surgical set should be used.

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If the products are not to be taken directly with hand instruments, they should be opened on the sterile cloth.

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Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacuum	Pre-Vacuum	3 times < 60 mbar
Temperature	134°C / /270°F	Minimum Drying time	10 minutes in dry area

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3.m Universal Straight Abutment: Universal Straight Abutment are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. It is milled and made suitable for the patient. It is fixed to the implant with the help of a screw. When preparing the universal abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant completely or the connection will not be provided. After the universal abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the universal abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. The milled universal abutments should be washed at 55 - 60 °C for 15 minutes. The products should not be sterilised and used without washing. Single use only.

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3.r Multi-Unit Tibase Abutment: Multi-Unit Tibase abutment; These are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants in completely edentulous patients. Multi-Unit Tibase abutment is designed as hex and no-hex. Both abutments are compatible with Venuscon Implants. Single use only.

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3.v Impression Parts: These are the parts used to determine the dimensions of the crown planned to be placed on the implant and superstructure. There are two separate options for classical and digital methods. Single use only.

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3.y Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool:

It is used to ensure that the Screw Retain Drill and Broken Screw Extractor Hand Tool, which are used to remove the broken screw in the mouth, can centre the screw to be removed. The Broken Screw Removal Centring Part is first placed in the area where the screw to be removed and then the screw is removed from the mouth with the Screw Retain Drill or Broken Screw Extractor Hand Tool. It is used for removing broken screws in the patient's mouth. After placement of the Broken Screw Removal Centering Part, the broken screw is removed with the Broken Screw Extractor Hand Tool. It is released as non-sterile. It can be used maximum 15 times after cleaning and sterilisation.

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3.ab Hand Tools : During the placement of implants in the mouth; It is used with carrying and clamping parts. It can be used up to 15 times after cleaning and sterilisation.

3.ac Gingival Height Part: The gingival height part is used to measure the level of the gingiva after the implant is placed into the bone. The gingival levels of the gingiva former and abutments to be used are determined after using the gingival height parts. It can be used after cleaning and sterilisation.

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Annex 8 – MKP.02-L.12 Medisolaris Impression Part-Screw Compatibility List

5. Sterilization and Packaging: Superstructure Components, Gingiva Formers, Prosthetic Parts and Instruments are packaged unsterilised, sealed in aluminium foil. All system parts must be sterilised before use on patients. All parts must be sterilised by autoclave according to the recommended protocols. The recommended steam sterilisation parameters are as in the table below.

Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacum	Pre-Vacum	3 times < 60 mbar
Temperature	134°C // 270°F	Minimum Drying Time	10 minutes indoor

The packaging of Superstructure Components, Gingiva Formers, Prosthetic Parts, Auxiliary Parts and Hand Instruments are not suitable for autoclave sterilisation. The products must be unpacked in a controlled environment before sterilisation, placed in an autoclave bag and then sterilisation must be performed.

6. Placement: Medisolaris dental abutments are placed on the respective implants with carrying and clamping parts and hex keys. Medisolaris ball abutments are self-screwing and are fixed to the implants in the hex slot above the round part using a ball

driver. Recommended torque values for screws and superstructures can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants.**

Cement should not be used in multi-unit abutment application. Ball abutments and Locator abutments should be used for removable prosthesis for patients who are used to removable prosthesis. Angled options of abutments should be used according to the posture of the implant.

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► Soft tissue health, ► Parafunctions. Venüscon Implants have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment.

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Not to use surgical hand instruments more than the maximum number of uses

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12. MEANING OF SYMBOLS

	Batch Code		Do not re-use		Sterilized using irradiation
	Katalog Numarası Catalogue Number		Do not re-sterilize		CE Marking / Notified Body Number
	Used by date		Date of manufacture		Non-Steril
	Do not use if package is damaged and consult instructions for use		Caution		UDI Number
	Consult instructions for use or consult electronic instructions for use		Temperature limit		Manufacturer
	Humidity imitation		Keep away from sunlight		Double sterile barrier System
	Medical Device		Unique Device Identifier		

MEDİSOLARIS SAĞLIK HİZMETLERİ SANAYİ ve TİCARET LİMİTED ŞİRKETİ - Fatih Mahallesi 1191. Sokak No: 28/1 B Blok
Apt. Gaziemir / İZMİR - +90 232 203 6655 / www.medisolaris.com

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www.medisolaris.com

A. Medisolaris Dental Implants

There are two separate series of Medisolaris Dental Implants, Merkurhex and Venuscon. Medisolaris Dental Implants are bone-level dental implants that are surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations.

A.1 Medisolaris Dental Implants Features

Medisolaris Dental Implant System is used for the treatment of mandibular and maxillary; single and partial edentulism.

A device specially designed to be surgically implanted in or over the mandibular or maxillary bone as a means of support and retention against dislodgement of the dental prosthesis. It is either threaded through the gum by the dentist (part of the implant emerges through the gum for direct support) or completely embedded under the gum (for removable denture support only).

Medisolaris Dental Implant is a bone-level dental implant that is surgically placed in the upper or lower jaw bone to support artificial dentures used to restore the patient's chewing function. From top to bottom, the implant's core geometry descends conically. Implants with a tapered core have higher tensile force and insertion torque. Reducing the core diameter without increasing the outer diameter too much increases the contact area between the bone implant and increases the tensile strength.

Diameter and Length Options

- Merkurhex Implant

Lenght \ Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.2	Available	Available	Available	Available	Available	Available
Ø 4.7	Available	Available	Available	Available	Available	Available
Ø 5.2	Available	Available	Available	Available	Available	Available

- Venuscon Implants

Lenght \ Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.3	Available	Available	Available	Available	Available	Available
Ø 5.0	Available	Available	Available	Available	Available	Available

- 1.a Warnings :** The information presented in the Instruction for Use is not sufficient for the immediate application of the Medisolaris Dental Implant System. The user must be legally qualified to perform surgical operations and must have advanced training in oral implantation and detailed knowledge of the system components before use. Surgical instruments must be sterilized (according to the manufacturer's recommendation) before use. Surgical instruments should not be used more than the maximum number of times, and the location of the implant socket should be adjusted correctly. There should be at least 1.5 mm between the intact teeth and at least 3 mm with the dental implant. The amount of bone in the vestibule, lingual or palatal should also be at least 1 mm. The bone temperature should not exceed 47C° during the preparation of the implant bed. Medisolaris Dental Implants have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. They have not been tested for heating, displacement or image artifact in the MR environment. The safety of Medisolaris Dental Implants in the MR environment is unknown. Scanning a patient with this device may result in patient injury.
- 1.b Legal Liability Disclosure:** The use of this product is the sole responsibility of the end user. No liability is assumed for loss or damage caused by misuse of this product, abnormal use (careless handling, sabotage and ignoring safety information), failure to comply with this IFU. This system can only be used with original Medisolaris Dental parts. The use of third party parts renders Medisolaris Sağlık Hizmetleri A.Ş. products ineffective.
- 1.c Copyright and Patent:** The logos of Medisolaris Dental Implants are protected patents and owned exclusively by Medisolaris Sağlık Hizmetleri A.Ş. any use of these marks is prohibited and is subject to prosecution under Turkish and international law.
- 1.d Education:** Dental implant treatment involves complex surgical and restorative procedures. Such procedures should be performed by dentists who have received the necessary training in this field. Inadequate training can lead to surgical and prosthetic complications, resulting in failure of the implant and restoration.
- 1.e Patient Populations:** Patients with an Implant Indication and suitable for the implant target patient population.
- 1.f Raw Material :** Stainless Steel
- 1.g Medisolaris Surgical Kit Usage Indications :** In implant-indicated patients, creating the implant socket, placing the implant in the socket and tightening the superstructures and screws.
- 1.h Contraindications :**
- Hard or soft tissue pathology at the implant site (untreated)
 - Active periodontal infection or untreated periodontitis
 - Local inflammation or mucosal diseases such as lichen planus
 - Patients with uncontrolled systemic problems that may negatively affect wound healing
 - Immunosuppressive disease
 - History of radiation therapy in the head or neck region
 - Used/using some intravenous (IV) bisphosphonates
 - Bone resorption and/or a disease condition affecting bone metabolism
 - History of alcohol and drug addiction
 - History of allergy to titanium
 - Have any impairment that interferes with the ability to maintain adequate daily oral hygiene
 - Patients with recent myocardial infarction and cerebrovascular events
- **Residual Risks :** Infection, injury during surgical procedure, perforation and parasthesia after surgical intervention, inability to regain chewing function or delay in regaining chewing function, cutting the gum and



hurting the patient, swallowing the implant, swallowing the Medisolaris drill extensions, swallowing the Medisolaris drill, swallowing the clamping piece, the need for a re-operation, Failure to create an implant socket, resulting in the product not serving the purpose of an artificial tooth root, creating stress in the bone - necrosis (dead) bone formation and failure to ensure osseointegration of the implant, prolongation of the operation, Periimplantitis (Implant), Perimucositis (Implant), Implant mobility (Implant), Failure to achieve osseointegration (Implant), Primary Stability (Implant), Inappropriate level of stable bone and soft tissue (Implant)

1.i Potential Side Effects:

• Infection	• Difficulty opening the mouth
• Bleeding and hematoma	• Tissue numbness
• Injuries to the nerves, sinus and nasal cavity	• Tingling in tissues
• Incorrect positioning of the implant	• Sensory Alteration
• Incision line patency	• Damage to neighboring tooth
• Bruising	• Sinus and nasal floor injury in the upper jaw
• Swelling	

1.j Placement Configuration of Implants in the Mouth:

Small diameter implants are placed in areas where the alveolar crest is thinner because the tooth itself and the root are small in diameter. After placing an implant in the alveolar crest, at least two mm of bone distance should remain around the implant in the bucco-lingual and mesio-distal directions. This two mm distance provides the necessary support to the implant during the use of the implant. In Venuscon and Merkurhex series implants; 3.2 mm diameter implants with a narrow platform are recommended to be placed in the area marked in blue on the side. It is recommended that the implant sizes to be used should be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental Implants can be used for all bone types. The implants placed in these areas are as follows. For recommended product codes, please refer to MKK.01 and MKK.03 IFU.

Likewise, large diameter implants are also used in areas where the diameter of the tooth and bone is thick, that is, in areas where there are teeth with more than one root. Large-diameter implants are placed in the area of the lower and upper jaw incisors where there is a biting force instead of chewing. In Medisolaris implants; $\varnothing 3.7/4.2/4.7/5.2$ mm diameter implants for the Merkurhex series with a regular platform and $\varnothing 3.7/4.3/5.0$ implants for the Venuscon series are recommended to be placed in the area marked in blue on the side. The implant sizes to be used are recommended to be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental Implants can be used for all bone types. The implants placed in these areas are as follows. For recommended product codes, please refer to MKK.01 and MKK.03 IFU.

1.k Bone Type

 <p>D4</p>	 <p>D3</p>	 <p>D2</p>	 <p>D1</p>
Thin trabecular bone	Thin trabecular bone with thinner porous cortical bone	Dense trabecular bone with porous cortical bone	Dense cortical bone

This surgical guide has been prepared for the various bone densities studied during the use of drills. The dentist is responsible for selecting the protocols to be applied according to the bone quality and anatomical structure during the study.

Bone quality is categorized by bone density. Bone quality is divided into 4 groups as D1, D2, D3 and D4. Bone quality is of great importance in implant dental treatment. It affects many factors such as implant planning, number of implants, surface properties of the implant, implant surgery technique, duration of implant treatment.

D1 bone is the bone quality group in which the bone is the densest and hardest. It is usually seen in the anterior region of the lower jaw, in cases where the bone has melted excessively. It is a group that should be considered in implant treatment. The bone is very hard. Care should be taken not to heat the bone during implant surgery. Otherwise, necrosis of the bone and loss of the implant is inevitable. D1 Bones may have less blood supply, therefore, implant - bone healing will not occur and implants may be lost.

D2 bone is the ideal bone for dental implants. The density and blood supply of the bone and the surface contact between the implant and bone are the best compared to the other groups.

The density of D3 bone is less dense than the density of D2 bone. It is the bone quality that should be approached more cautiously when applying implants.

D4 bone is the bone quality where the density of the bone is close to zero.

In D4 bone quality, a stable fit of the dental implant to the jawbone may not be achieved. In such a case, it is recommended to place the implants in the area to be implanted, after opening the implant socket, compressing implant bone powder into the socket and waiting for 3-6 months.

1.l Opening of Steril Dental Implant Package

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STEP 2: Open the back of the blister. The product should not be used if the adhesion area is torn, broken or opened. The implant should be opened in a sterile environment where the surgical operation will be performed.

STEP 3: If the implant is to be taken with anguldurva; the key suitable for the anguldurva is attached to the tip of the anguldurva. If the implant is to be taken manually, the manual implant key is inserted into the torque tool.

STEP 4: The cap of the tube is held upwards. After the cap of the tube is turned open, the implant in the stand is removed with the anguldurvaimplant key or manual implant key.

STEP 5: To remove the implant from the stand, the anguldurvaimplant key or manual implant key is inserted into the implant. Make sure that the keys are fully inserted into the implant so that the implant does not fall out.

STEP 6: When the implants are removed from the carrier, they are removed by moving back and forth.

STEP 7: If the anguldurva is used to place the implant, the appropriate torque is adjusted from the motor and placed in the slot opened in the bone. If the implant is placed manually, the implant is removed from the socket and placed in the opened socket after the appropriate torque is adjusted before taking the implant with the torque.

STEP 8: To place a cap on the implant, the cap in the tube is removed in the surgical operation room with the help of a key suitable for the cap.

STEP 9: The cap removed from the tube is placed on the implant without manual contact. If the products are not to be taken directly with hand tools, they should be opened on the sterile cloth.

B. General Characteristics of Drills

All drills are compatible with Venuscon and Merkurhex implants. Drills should be selected according to the diameter and length of the implant to be placed. Among the drills designed for Venuscon and Merkurhex implant system, pilot drills are designed for every size implant with a stopper. With the help of the length lines on the final drills, it is possible to open the appropriate slot for the implant. When the drills are used without a stopper, there are black guide lines on the drills to determine the slot depth. Each of these lines represents a specific length. The lengths they represent are shown in the picture below. The dentist should open the slot using the guide lines according to the length of the implant to be placed. In drills designed for Merkurhex and Venuscon implant system, it is appropriate to use the next drill for hard bones.

All surgical instruments are made of stainless steel.

It is recommended to use each driller maximum 15 times. Before each use, the drills should be cleaned in accordance with the instructions under A.1 General Cleaning Recommendations and sterilized in accordance with the parameters under A.2 Sterilization Process.

Drills should be stored in accordance with A.3 Storage conditions.

B.1 General Cleaning Recommendations

- Observe general precautions for handling contaminated and biohazardous materials.
- Clean immediately after each use to prevent biological fluids and tissues from drying on the instruments.
- Whenever possible, disassemble parts and tools before cleaning.

- Do not rely solely on automatic cleaning. Detailed manual cleaning is recommended.
- Pre-cleaning should include wiping parts, soaking in warm enzymatic solution for at least twenty (20) minutes, and rinsing with tap water.
- Routine cleaning should include (a) washing parts using a broad spectrum cleaning solution, followed by rinsing and
- (b) parts fully immersed in the cleaning solution for at least ten (10) minutes, preferably
- It should include selenizing at 45-50kHz followed by thorough rinsing and drying.
- Dry immediately and completely to avoid oxidation.

B.2 Sterilization Process

Surgical Hand Tools must be sterilized by autoclave method before use. In autoclave sterilization (steam), they should be sterilized at 134°C for 3 minutes.

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The packaging of Surgical Hand Tools is not suitable for autoclave sterilization. The products must be unpacked in a controlled environment before sterilization, placed in a package suitable for sterilization and then the sterilization process must be performed.

* Make sure that the autoclave device to be sterilized is maintained annually and its parameters are correct.

B.3 Storage Conditions

It is recommended that the areas where Surgical Hand Instruments are stored comply with the following temperatures and humidity values.

Temperature	18 C° - 28 C°
Humidity	%40 - %60

B.4 Surgical Instruments

- 1. Punta Drill :** The punta drill is the initial drill and is used to determine where the implant will be placed. There is no coating on the surface. A heat treatment process is applied. The Punto drill is supplied non-sterile and must be sterilized by the end user.
- 2. Round Drill :** It is the drill used after the punto drill and is used for widening and cleaning the implant placement site after the point. There is no coating on the surface. Heat treatment process is applied. Round drill is supplied non-sterile and must be sterilized by the end user.
- 3. Merkurhex ve Venuscon Pilot and Final Drills :** The drills are osteotomy drills. Osteotomy drills are used to make a suitable slot in the bone for implant placement. Pilot drills are designed with a stopper to increase safety of use. Pilot drills are 6/8/10/11.5/13/16 mm in length. Pilot drill will be used in soft bones for implants with Ø3.2 diameter. Merkurhex and Venuscon Final Drills are available in five different diameters and 2 lengths. Final drills are designed in 2 types: short and long. Short drills are used where the mouth opening is small and long drills are used where the mouth opening is sufficient. Short drills are used to open slots from 6 mm to 10 mm. Long drills are suitable for slotting up to 16 mm. Drills are designed with two wings. Codes are written on them with laser marking according to the implant diameter to be used. All drills are designed for motorized use. There is no coating on the surface. Heat treatment process is applied. Sharpening process is applied. Drills are supplied non-sterile and must be sterilized by the end user.

4. **Countersink Drill** : In cases where the patient's bone is too hard, neck expander drills are used after the final drills to expand the cortical bone and reduce the stress on the neck area of the implant. Countersink drills are available in five different diameters and one size. The countersink drills have coloring and guide lines to guide the diameter and length of the drills. All neck expander drills are designed for motorized use. There is no coating on the surface. Heat treatment process has been applied. Countersink drills are supplied non-sterile and must be sterilized by the end user.
5. **Drill Extension** : Used for extending drill lengths. There is no coating on the surface. Heat treatment process is applied. Drill extension is supplied non-sterile and must be sterilized by the end user.
6. **Rachet Adaptor** : It is an apparatus that allows the use of surgical parts used with a motor. There is no coating on the surface. Heat treatment process is applied. The rchette adaptor is supplied non-sterile and must be sterilized by the end user.
7. **Transferring and Tightening Part** : Implant Transferring and Tightening Parts are used for placing implants in the implant socket. Short clamping and handling parts are used where the mouth opening is small; long clamping and handling parts are used where the mouth opening is larger.
8. **Screw Retain Drill**: Screw retain drill are drills used to remove broken implants from the socket in the bone.

B.5 Preparation of the implant socket

In order to place Medisolaris Dental Implants in the bone, the implant socket should be prepared by following the sequence described in the image below. This protocol applies to soft bones and hard bones. Medisolaris final drills have a single length. In order to open the slot according to the implant size, it is recommended to open the slot by following the depth lines shown by marking after using pilot drills with special stoppers for each size.

To prevent bone heating and burning, it is recommended to use the drills with a physio dispenser that can be supplemented with saline fluid. It should be used after making sure that the drills are fully seated in the angul durwa connected to the physio dispenser. It is recommended to apply with 700 rpm.

It is recommended to use each driller a maximum of 15 times. Each implant should be considered as one use.

A.5.a Recommended Torq Values

Process	Recommended Torq Values
Placement of Implants	35-40 Ncm
Placement of Gingiva Formers	15 Ncm
Placement of Covers	15 Ncm
Tightening of Merkurhex Screw	30 Ncm
Tightening of Venuscon Screw	30 Ncm
Tightening of Multi-Unit Oclusal Screw	20 Ncm
Tightening of Screw Abutments	30 Ncm

The protocol to be followed for the Venuscon implant is as follows.

Venuscon | 3.2mm

3.0/3.2

DRL 2023-16
2.0/2.3mm

Final Drill

DRL 2328
2.3/2.8mm

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CBD 3200

Venuscon | 3.7mm

3.0/3.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

Final Drill

DRL 2834
2.8/3.4mm

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CBD 3700

Venuscon | 4.3mm

3.0/4.3

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

Final Drill

DRL 3438
3.4/3.8mm

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CBD 4200

Venuscon | 5.0mm

3.0/5.0

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

Final Drill

DRL 4146
4.1/4.6mm

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CBD 5000

The protocol to be followed for the Merkurhex implant is as follows.

Merkürhex | 3.2mm

3.0/3.2

DRL 2023-16
2.0/2.3mm

Final Drill

DRL 2328
2.3/2.8mm

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CBD 3200

Merkürhex | 3.7mm

3.5/3.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

Final Drill

DRL 2834
2.8/3.4mm

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CBD 3700

Merkürhex | 4.2mm

3.5/4.2

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

Final Drill

DRL 3438
3.4/3.8mm

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CBD 4200

Merkürhex | 4.7mm

3.5/4.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

Final Drill

DRL 3844
3.8/4.4mm

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CBD 4700

Merkürhex | 5.2mm

3.5/5.2

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

DRL 3844
3.8/4.4mm

Final Drill

DRL 4448
4.4/4.8mm

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CBD 5200

B.5 Clinical Benefit : Dental implants and Surgical Instruments are a safe, functional alternative treatment. For example, it significantly improves the patient's quality of life, improves the patient's chewing function, greatly reduces loss of function and provides symptom relief. It is in the form of natural teeth, provides chewing function, a better aesthetic appearance, healthy nutrition and self-confidence with healthy teeth.

A strong and stable dental implant restores the missing tooth so that it looks and functions like a natural tooth. Dental implants are advantageous because they can improve the stability of removable dentures or address the need for fixed or removable dentures. In a study by Perrotti et al, 2 different implant systems were found to be suitable for long-term treatment, with stable

marginal bone levels and a 10-year success rate. Dental implants help you maintain the natural shape of your face and smile. Leaving a gap in your mouth after losing one or more teeth leads to further health problems such as jawbone deterioration and loss. When the natural tooth is not supported, the jawbone deteriorates and loses its strength and stability. Dental implants are the only dental restoration option that preserves and stimulates natural bone, helps bone to grow and prevents bone loss.

B.6 Safety and Performance Information: Test reports that meet the safety and performance criteria of Medisolaris products have been applied and the results of the test reports have been found to be appropriate. Negative incidents that may occur in Medisolaris products should be reported to the relevant competent authority.

B.7 Disposal: Implants removed from the patient should not be re-implanted as their fatigue condition cannot be visually detected. These implants should be considered as biocontaminated and should be treated accordingly. Removed implants should be treated according to their biological hazards and disposed of in accordance with the legal regulations of the country and hospital policies.

B.8 Meaning of Symbols :

	Batch code		Do not re-use		Sterilized using irradiation
	Catalogue number		Do not re-sterilize		CE Marking / Notified Body number
	Used by date		Date of manufacture		Non-Steril
	Do not use if package is damaged and consult instructions for use		Caution		UDI Number
	Consult instructions for use or consult electronic instructions for use		Temperature limit		Manufacturer
	Humidity imitation		Keep away from sunlight		Double sterile barrier system
	Medical Device		Unique Device Identifier		

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1-) DENTAL IMPLANTS

Medisolaris (MERKURHEX) Dental Implant and Superstructures produced by MEDİSOLARİS Sağlık Hizmetleri San. Tic. Ltd. It is and is surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations. MERKURHEX Dental Implants are bone level dental implants. The target users for Venuscon Dental Implants are dentists.

Superstructures are the device that used to support crowns, bridges or superstructures in the lower and/or upper jaw, used with Dental Implants in completely edentulous or partially edentulous patients. It is recommended that Medisolaris products not be used for purposes other than their stated intended use.

Medisolaris Dental Implant Basic Information :
For Merkurhex Series Implant :

Lenght : 6 / 8 / 10 / 11.5 / 13 / 16 mm

Diameter : 3.2/3.7/4.2/4.7/5.2 mm

Lenght Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.2	Available	Available	Available	Available	Available	Available
Ø 4.7	Available	Available	Available	Available	Available	Available
Ø 5.2	Available	Available	Available	Available	Available	Available

1.a Warning : The information presented here is not sufficient for the immediate application of the Medisolaris Dental Implant System. The user must be legally competent to perform surgical operations and must receive advanced training in the field of oral implantation and detailed knowledge of the system components before use.

1.b Legal Liability Disclosure: The use of this product is the sole responsibility of the end user. No liability is assumed for loss or damage caused by misuse of this product, abnormal use (careless handling, sabotage and ignoring safety information), failure to comply with this IFU. This system can only be used with original Medisolaris Dental parts. The use of third party parts renders Medisolaris Sağlık Hizmetleri A.Ş. products ineffective.

1.c Copyright and Patent: The logos of Medisolaris Dental Implants are protected patents and owned exclusively by Medisolaris Sağlık Hizmetleri A.Ş. any use of these marks is prohibited and is subject to prosecution under Turkish and international law.

1.d Education: Dental implant treatment involves complex surgical and restorative procedures. Such procedures should be performed by dentists who have received the necessary training in this field. Inadequate training can lead to surgical and prosthetic complications, resulting in failure of the implant and restoration.

1.e Patient Population : Patients with complete bone development (18 years of age or older), not pregnant, not in the rapid puberty period of craniofacial changes and not during lactation, without hard or soft tissue pathology in the implant area, without active periodontal infection or untreated periodontitis, after periodontal treatment; Good oral hygiene, good response to periodontal treatment, non/smoker, systemically healthy patients, no unstable endocrine diseases (such as diabetes), no local inflammation or mucosal diseases such as lichen planus, no osteoporosis (Bone resorption), no immunosuppressive disease, no cancer cases requiring chemotherapy / radiotherapy / anti-resorptive drug use, no acquired or acquired coagulation disorders; no intravenous (IV) based bisphosphonate use / use, no history of alcohol and drug addiction, no history of allergy to titanium and stainless steel, no systemic disease requiring regular medication, no mental development disorders such as mental retardation. Venuscon Implants have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. They have not been tested for heating, displacement or image artefact in the MR environment. Venuscon safety in the MR environment is unknown. Scanning a patient with this device may result in patient injury.

1.f Raw Materials and Ingredients: Implants and superstructures are made of ISO 5832-3 Ti6Al4V ELI / GRADE23, surgical instruments are made of 1.4116, hand instruments are made of 1.4057 Stainless Steel. The surface of Medisolaris Dental Implant is roughened by sandblasting and passivation SBM method consisting of hydroxyapatite and tricalcium phosphate to obtain surface microstructure, osseointegration, Ra values 2-2.5Ra.

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- Total and partially edentulous patients
- Mandible and maxillary defects
- In patients who have difficulties in the use of removable partial dentures
- In the presence of insufficient number and position of support teeth for removable and fixed prosthesis
- In order to prevent the cutting of the teeth adjacent to the missing tooth, to support the bridge prosthesis used in the treatment of missing teeth.

1.h Contraindications:

✓ Having hard or soft tissue pathology in the Implant area	✓ Those with active periodontal infection or untreated periodontitis
✓ Those with local inflammation or mucosal diseases such as lichen planus	✓ Patients with systemic problems that cannot be controlled and will negatively affect wound healing
✓ People with immunosuppressive disease or on medication	✓ A disease that affects osteoporosis and/or bone metabolism
✓ Have a history of radiation therapy to the head or neck area	✓ Having a history of alcohol and drug addiction
✓ Used/used intravenous (IV) based bisphosphonate	✓ Have a history of allergy to titanium
✓ Patients with recent myocardial infarction and cerebrovascular accident	✓ Have any disorder that interferes with the ability to maintain adequate daily oral hygiene or plaque control

1.i Residual Risks:


• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening

1.j Potential Side Effects:

• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism

• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening


1.k Placement Configuration of Implants in the Mouth: Small diameter implants are placed in areas where the alveolar crest is

thinner due to the small diameter of the tooth itself and the root. After the implant is placed in the alveolar crest, at least two mm bone distance should remain around the implant in the bucco-lingual and mesio-distal directions. This two mm distance provides the necessary support to the implant during the use of the implant. In Merckürhex series implants; 3.2 mm diameter implants with a narrow platform are recommended to be placed in the area marked in dark colour on the side. It is recommended that the implant lengths to be used should be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone

types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME
SLM 3208	Merkürhex Implant Ø3.2mm L:8mm Platform 3.0mm
SLM 3210	Merkürhex Implant Ø3.2 mm L:10 mm Platform 3.0 mm
SLM 3211	Merkürhex Implant Ø3.2 mm L:11.5 mm Platform 3.0 mm
SLM 3213	Merkürhex Implant Ø3.2 mm L:13 mm Platform 3.0 mm
SLM 3216	Merkürhex Implant Ø3.2 mm L:16 mm Platform 3.0 mm

Likewise, large diameter implants are also used in areas where the diameter of the tooth and bone is thick, that is, where there are teeth with more than one root. A large diameter implant is placed in the area where the biting force is located instead of chewing, where the lower and upper jaw incisors are located. In Medisolaris implants, it is recommended to place implants with Ø3.7/4.2/4.7/5.2 mm diameter in the area marked with dark colour on the side for Merckürhex series with regular platform. The implant sizes to be used are recommended to be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME	PRODUCT CODE	PRODUCT NAME
SLM 3706	Merkürhex Implant Ø3.7mm L:6mm Platform 3.5mm	SLM 4206	Merkürhex Implant Ø4.2mm L:6mm Platform 3.5mm
SLM 3708	Merkürhex Implant Ø3.7mm L:8mm Platform 3.5 mm	SLM 4208	Merkürhex Implant Ø4.2mm L:8mm Platform 3.5mm
SLM 3710	Merkürhex Implant Ø3.7mm L:10mm Platform 3.5 mm	SLM 4210	Merkürhex Implant Ø4.2mm L:10mm Platform 3.5mm
SLM 3711	Merkürhex Implant Ø3.7mm L:11.5mm Platform 3.5mm	SLM 4211	Merkürhex Implant Ø4.2mm L:11.5mm Platform 3.5mm
SLM 3713	Merkürhex Implant Ø3.7 mm L:13 mm Platform 3.5 mm	SLM 4213	Merkürhex Implant Ø4.2mm L:1 mm Platform 3.5mm
SLM 3716	Merkürhex Implant Ø3.7 mm L:16 mm Platform 3.5 mm	SLM 4216	Merkürhex Implant Ø4.2mm L:16mm Platform 3.5mm
SLM 4706	Merkürhex Implant Ø4.7 mm L:6 mm Platform 3.5 mm	SLM 5206	Merkürhex Implant Ø5.2mm L:6mm Platform 3.5mm
SLM 4708	Merkürhex Implant Ø4.7 mm L:8 mm Platform 3.5 mm	SLM 5208	Merkürhex Implant Ø5.2mm L:8mm Platform 3.5mm
SLM 4710	Merkürhex Implant Ø4.7 mm L:10 mm Platform 3.5 mm	SLM 5210	Merkürhex Implant Ø5.2 mm L:10 mm Platform 3.5 mm
SLM 4711	Merkürhex Implant Ø4.7 mm L:11.5 mm Platform 3.5mm	SLM 5211	Merkürhex Implant Ø5.2 mm L:11.5 mm Platform 3.5 mm
SLM 4713	Merkürhex Implant Ø4.7 mm L:13 mm Platform 3.5 mm	SLM 5213	Merkürhex Implant Ø5.2 mm L:13 mm Platform 3.5 mm
SLM 4716	Merkürhex Implant Ø4.7 mm L:16 mm Platform 3.5 mm	SLM 5216	Merkürhex Implant Ø5.2 mm L:16 mm Platform 3.5 mm

1.l Packaging and Sterilization : Medisolaris Dental implants are packaged sterile with Caps. The products are disposable. Gamma radiation is used to sterilize implants after packaging. Thanks to sealed packaging, implants are prevented from being affected by external environments. For proper storage conditions, the temperature should be between 18-28°C and the humidity should be between 40-60%. Sterilization is valid until the expiry time stated on the label on the package. Implants in damaged sterile packages should not be used. Contaminated implants cannot be resterilized. Reuse of implants may result in clinical failure and complications including, but not limited to, infection, failure to osseointegrate implants, bone loss, and damage to adjacent structures. The implant surface should not be touched by hand, as this may deteriorate its surface properties. All other superstructures and hand tools should be sterilized by autoclave. For recommended autoclave sterilization parameters, see the Superstructure section and *MKK.02 Surgical Instruction for Use*. Do not use the product whose packaging is damaged.

1.m Shelf Life: The expiration date is clearly stated on each implant package.

1.n Treatment Planning: Comprehensive patient examination, pre-operative imaging, pre-operative diagnosis and treatment planning should be done before implant treatment. In patients who cannot achieve the required bone density for the implant, it is recommended to apply for bone augmentation procedures with advanced surgical techniques to increase bone density before the application and then perform the application. Inadequate planning and evaluation can lead to implant complications or even implant failure.

1.o Preparation of the Surgical Area: Preparation of the implant area should be done with the surgical drills in the Medisolaris Dental surgical set. The drilling speed should be decided by the physician according to the hardness level of the bone. Drilling should be performed under abundant irrigation. Care should be taken to prevent the bone from getting too hot. The use of surgical drills that have lost their sharpness should be avoided. Because these may cause unnecessary trauma to the bone. It is recommended to follow the drilling sequence according to the implant diameter to be used. For detailed information about the surgical operation, see *MKK.02 Surgical Instruction for Use*.

1.p Surgical Placement of Implants: Check the length and diameter of the implant on the product label. Attach the (3) labels included in the packaging to the patient's file. Open the outer cardboard package and remove the thermoformed packaging. Open the tyvek paper on the thermoform and empty the sterile tube into the surgical area. The implant tube consists of a single part: Open the implant tube by turning the top cap. Use the implant carrier in the set to move the implant into the surgically prepared implant slot. The ratchets or physiodispenser included in the set should be used to place the implant. Medisolaris Covers should be closed with the help of a 1.25 mm hex key. Caps and healing caps should be torqued by hand.

1.q Recommended Torq Values

Process	Recommended Torq Values	Process	Recommended Torq Values
Placement of Implants	35-40 Ncm	Screw Venuscon	30 Ncm
Gingiva Former	15 Ncm	Multi-Unit occlusal	20 Ncm
Cover	15 Ncm	Self-Screwing Abutments	30 Ncm

For detailed information see *MKK.02 Instruction for Use*. *Recommended used at 700 RPM.

1.r Healing Process of Medisolaris Dental Implants: In general, implants should be left to heal for 2 to 4 months, depending on the patient's bone quality, type and general health status. Immediate loading can be performed with the right case selection and adherence to the immediate loading protocols accepted in the literature. If immediate loading is to be performed, the occlusion environment should be checked. If there are multiple implants, they should be splinted and in cases of complete edentulism, at least 4 implants in the mandible and 6 implants in the maxilla should be connected to each other with a prosthesis. Depending on the bone anatomy, the appropriate number, diameter and length of implants should be placed.

1.s Information to be communicated to the patient after implant placement

It is important that patients follow the instructions given below after dental implant surgery to protect their implants for life. For detailed information, please refer to *MKK.04 Patient Education Brochure*. It is recommended to be given to patients after the operation.

In dental implant surgical treatment, it is recommended to wait at least 4 months in the upper jaw and at least 3 months in the lower jaw for the implant to fuse with the jaw bone.

After the dental implant is placed, the patient will feel some pain and discomfort. The pain level is usually mild. After the operation, the dentist may prescribe painkillers. With painkillers, the pain can be controlled. These should be taken regularly for 2 days. The first dose of painkillers must be taken before the local anaesthetic wears off. The pain will begin to subside within a few days. If the pain does not go away after a few days, you should consult your doctor.

1. After the dental implant is placed, the patient will feel some pain and discomfort. The pain level is usually mild. After the surgery, the dentist may prescribe painkillers. With painkillers, the pain can be controlled. These should be taken regularly for 2 days. The first dose of painkillers should be taken before the local anesthetic wears off. The pain will start to decrease within a few days. If the pain does not go away after a few days, you should consult your doctor.

2. After the dental implant, the dentist may prescribe antibiotics and antiseptic mouthwash to the patient if deemed appropriate.
3. A sterile sponge is placed on the area where dental implant surgery is performed to control bleeding. The sterile sponge should be bitten by the patient with light force for 20 minutes. At the end of 20 minutes, the sponge should be removed.
4. The patient should not eat or drink anything for two hours after the operation. The patient should not spit or rinse his/her mouth with any liquid for two hours. After two hours, the patient should not eat hard and hot foods on the first day. Hot foods increase both bleeding and the risk of infection.
5. It is normal to have bleeding in the form of leakage during the first 24 hours after the procedure. To control bleeding, the patient bites gauze for 30 minutes. If the bleeding continues for more than 2 days, you should consult our physician. Aspirin and aspirin compounds should not be taken as painkillers because they cause bleeding tendency.
6. Swelling may occur after implant surgery. This is a normal condition. If the pain and swelling do not disappear after 5-6 days, you should consult your physician.
7. The patient should avoid tampering with the area where the dental implant is placed with the tongue, finger or hard instruments.
8. The patient should rest a little after the implant surgery and avoid strenuous exercise for 48 hours.
9. If the patient is using removable dentures, he/she should not use them until the date that the dentist deems appropriate.
10. In dental implant surgery, the wound surfaces are closed with sutures. The sutures must be removed by the dentist after one week at the earliest. The patient is informed by the dentist about the early period in which the treated tooth part can be used.
11. After dental implant surgery, the patient should continue oral care after the first 24 hours. The patient should use antiseptic mouthwash prescribed by the dentist three times a day for the first week after the surgery. Antiseptic mouthwash should not be used after one week. After antiseptic mouthwash, the patient can gargle with warm salt water. Add 1 teaspoon of salt to 1 glass of water and gargle 3-4 times a day to disinfect the area. This will help clean the healing area.
12. After 48 hours, the surgical site should be lightly brushed. Other teeth can be brushed and cleaned normally after the operation. There should be no food residue in the surgical site and the patient should rinse the mouth with warm water after meals.
13. Smoking delays the healing process and prevents dental implants from staying in the mouth for a long time, so the patient should stop smoking cigarettes and similar tobacco before starting dental implant treatment.
14. It is very important to follow all these instructions mentioned above for the success of dental implant surgery. If there are signs of infection after implant surgery or if the patient experiences nausea, fever or vomiting, he/she should consult his/her dentist.

1.t Long-term care of dental implants: In order to be able to use implants for a lifetime without any problems, it is important that they are properly cared for as described above.

1.u The importance of oral hygiene in dental implant patients: Implants are more susceptible to plaque and calculus build-up, inflammation and bone loss than natural teeth. In implant patients, good oral hygiene is very important. The patient should brush teeth correctly and clean the interdental and interimplant areas correctly. Oral hygiene kits are available for the private use of implant patients, which can be used by the patient. If the patient follows the oral hygiene instructions, inflammation of the periimplant tissues will be minimal. Removal of plaque and calculus deposits from the teeth will reduce the risk of periimplantitis, the most common cause of implant failure.

1.v Regular dental visits and examinations for implant patients: In patients treated with dental implant prostheses, check-ups are usually recommended every 3 months for the first 1-2 years and then every 6 months thereafter. At each visit, oral hygiene instructions are reinforced and the patient's ability to maintain good oral hygiene is checked. Patients with implant-retained removable prostheses are given some special instructions. The patient is told how to insert and remove the prosthesis without causing any damage to the implant or prosthesis. Excessive pressure should not be applied to the restoration placed on the implant. Restorations on implants are not as strong as natural teeth. The patient should avoid eating hard foods that can break the restoration. Avoid irritation of the implant area.

1.w Lifespan of dental implants: The lifespan of dental implants depends on many factors, such as patient health and the care provided by the patient. The lifespan of implants depends on how old the patient had the implants, how they live and how they care for them. For the success of the implants, the patient should maintain good oral hygiene and visit the dentist regularly for check-ups. Although a clear lifetime cannot be defined for the device, it is thought to have a lifetime of 10 years and longer under normal conditions. Implants have a shelf life of 5 years and superstructures do not have a shelf life. Implants are an effective solution for missing teeth and can be used for a lifetime with good care. Brushing, flossing and regular visits to the dentist after the implant prosthesis are important

to prolong the implant's stay in the mouth. Implants can replace a single tooth, a group of teeth or all upper and lower jaw teeth.

1.x Opening Sterile Packages: To take Medisolaris Dental implants from the product packaging with an angldurva, implant keys suitable for angldurva are used. To take the implants manually, the implant key suitable for manual removal is used. When taking the implants, the size suitable for the area to be placed in the mouth should be selected. For all diameters of the Venuscon implant, the 2.3 mm hexagon crimping spanner included in the surgical set should be used.

In order to attach the closing screws of Medisolaris Dental implants to the implants, the caps that are compatible in diameter that come out of the product are placed with a hex key.

STEP 1: Open the implant box and remove the sterile packaged thermoform.

STEP 2: The back of the thermoform is opened. If the adhesion area is torn, broken or opened, the product should not be used. The implant should be opened in a sterile environment where the surgical operation will be performed.

STEP 3: If the implant is to be taken with an angldurva; the key suitable for the angldurva is attached to the tip of the angldurva. If the implant is to be taken manually, the manual implant key is attached to the torque racket.

STEP 4: The cap of the tube is held upwards. After the cap of the tube is opened by turning the cap, the implant in the sling is removed with an angldurva implant key or manual implant key.

STEP 5: To remove the implant from the sling, the anguldruva implant key or manual implant key is inserted into the implant. Make sure that the keys are fully inserted into the implant so that the implant does not fall out.

STEP 6: When the implants are removed from the discs in the sling, they are removed by moving back and forth.

STEP 7: If an angldurva is used to place the implant, the appropriate torque is adjusted from the motor and placed in the slot opened in the bone. If the implant is placed manually, the implant is removed from the slot and placed in the opened slot after the appropriate torque is adjusted before taking the implant to the torque raschette.

STEP 8: In order to place a cover on the implant, the cover in the tube is removed with the help of a key suitable for the cover in the surgical operation room.

STEP 9: The cap removed from the tube is placed on the implant without manual contact.

If the products are not to be taken directly with hand tools, they should be opened on the sterile cloth.

1.y Post-Operative Care: Close follow-up of the patient after surgery is very important. It is recommended to use appropriate oral rinses and to ensure good oral care and not to use tobacco products.

2-) Medisolaris Dental Cover

The closure screw is used to cover and protect the implant after placing the implant in the bone. The closure screw can remain in the mouth for 2-3 months during the surgical healing period.

2.a Raw Material : ISO 5832-23 Ti6AL4V ELI

2.b Indications: The caps come out of the implant package used. It is used in two-stage implant surgery to prevent soft or hard tissue growth on the hollow inner surface of the implant.

2.c Contraindications: Contraindications are mentioned above.

2.d Possible Side Effects: Possible side effects are mentioned above.

2.e Placement: The caps are torqued by hand using a 1.25 mm hex spanner. Recommended torque values can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants.**

2.f Packaging and Sterilisation: The caps are packed together with the implant in a sterile tube. If the package is intact, sterilisation is valid until the expiry date on the package. In cases where the closure is single packaged, the product is single packaged without sterility and packaged with aluminium foil without sterility and must be sterilised by autoclave method before use. In autoclave sterilisation (steam), it should be sterilised at 134°C for 3 minutes.

Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacum	Pre-Vacuum	3 times < 60 mbar
Temperature	134°C / /270°F	Minimum Drying time	10 minutes in dry area

The packaging of the lids is not suitable for autoclave sterilisation. The products must be unpacked in a controlled environment before sterilisation, placed in a package suitable for sterilisation and then sterilisation must be performed.

3-) Medisolaris Dental Superstructures, Gingiva Formers, Prothetic Parts and Instruments

Superstructures are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with dental implants in completely edentulous or partially edentulous patients.

3.a Raw Material : ISO 5832-23 Ti6AL4V ELI- Stainless Steel 1.4116 ve 1.4057

3.b Gingiva Formers: Healing caps are used after the healing of the implant in the jawbone is completed and the overlying soft tissue is removed, until prosthetic treatment for gum shaping. Healing caps remain in the mouth for 15 days or less. They are for single use only.

3.c Temporary Abutment: Temporary abutment is used with Medisolaris implant in completely edentulous or partially edentulous patients. In immediate loading treatments, they are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw for more than 24 hours. Temporary abutment is an endosseous dental implant abutment designed for single use as a temporary prosthesis during the healing process while the permanent prosthesis is being produced. Temporary abutments are compatible with all diameters of Medisolaris implant for Venuscon series. It is fixed to the implant with the help of a screw.

3.d Straight Abutment: In order to use a cementable abutment, the implant must be placed in the patient without an angle. Cementable Abutments are available for both brands. It is compatible with all diameters of Medisolaris Dental implants for Venuscon series. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is bonded on it with the help of cement. Single use only.

3.e Angled Abutment: Angled abutments facilitate prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid anatomical structure when placing the implant. The angled abutment is compatible with all diameters of the Medisolaris implant for the Venuscon series. It is fixed to the implant with the help of a screw. The angled abutment is used in cemented restorations. It is for single use only.

3.f Anatomical Straight Abutment: The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This offers a more aesthetic option to the physician. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it emerges through the gum. In a possible gum recession, the underlying titanium is prevented from being seen. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zirconia is glued on it. It is disposable.

3.g Anatomically Angled Abutment : The use of an anatomically angled abutment facilitates prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid the anatomical structure when placing the implant. In addition, the use of aesthetic angled abutment shortens the healing time. The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This gives the dentist a more aesthetic option. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it comes out of the gum. In a possible gum recession, the underlying titanium is prevented from being seen. For each brand (Medisolaris Dental) implant, the same brand Anatomically Angled Abutment should be used. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is glued on it. It is for single use only.

3.h Multi Unit Abument: It is the external connection of the implant embedded flat into the bone. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different gum options are available. It is self-screwed and fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. It is disposable.

3.i Multi Unit Angled Abutment: It is the external connection of the implant embedded into the bone at an angle. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different angle and gum options are available. It is fixed to the implant with a screw. A tooth is made of materials such as porcelain or zircon. There are two types of multi-unit angled abutment compatible with Medisolaris implants: hex and non-hex. The hex structure is recommended to be used in single tooth restorations, while the non-hex structure is recommended to be used in bridge, bar and hybrid treatments. Hex structure is designed to prevent rotation in single tooth applications. Single use only.

3.j Ball Abutment: It is the external connection of the implant embedded flat into the bone. It is used in cases where the physician wants to apply removable prosthesis to the patient. Different diameters and gum lengths are available. It is self-screwed and fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. Single use only.

3.k Tibase Cad-Cam Abutment: It is the platform of the abutment inside the implant. Zirconia or custom (personalised) abutment is produced on it. It is produced for the purpose of using a company-specific platform. It is for single use only.

3.l Locater Abutment: It is an abutment system with low vertical height used for implant over implant removable prosthesis in completely edentulous jaws. The locater abutment piece is fixed to the implant with a metal cap and rubber in the removable prosthesis. Single use only.

3.m Universal Straight Abutment: Universal Straight Abutment are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. It is milled and made suitable for the patient. It is fixed to the implant with the help of a screw. When preparing the universal abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant completely or the connection will not be provided. After the universal abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the universal abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. The milled

universal abutments should be washed at 55 - 60 °C for 15 minutes. The products should not be sterilised and used without washing. Single use only.

3.n Premill Abutment: It is produced for CAD-CAM devices. The laboratory connects this superstructure to the device and makes a patient-specific abutment. While preparing the premill abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant properly or the connection will not be provided. When preparing the premill abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. Single use only.

3.o Ti-base Casting Abutment: For use in personalised cast prosthetic restorations in single or multiple implants. It is fixed to the implant with the help of a screw. Implants for Venuscon series in Medisolaris Dental implants are designed to be compatible with all superstructures. Interface connection type is the same. Single use only.

3.p Screws: The screw allows the superstructures to be mounted to the implant. The screws used with Measuring Parts and Analogues are different. It is packaged and sold together with the related product. Single use only.

3.q Titanium Cap: In completely or partially edentulous patients; These are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants. Titanium Caps are designed as Hex and No-Hex. Both caps are compatible with Venuscon Multi Unit Abutments. Ø 5.0 prosthetic platform diameter and 12.0 mm gingival height. No coating process is applied on the surface. Single use only.

3.r Multi-Unit Tibase Abutment : Multi-Unit Tibase abutment; These are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants in completely edentulous patients. Multi-Unit Tibase abutment is designed as hex and no-hex. Both abutments are compatible with Venuscon Implants. Single use only.

3.s Comfort Cap: The comfort cap prevents the patient from being uncomfortable by covering the multi unit abutment in the patient's mouth until the crown is loaded in Multi unit abutments. Comfort Caps are designed with 5.0 prosthetic platform diameter and 4.8 mm gingival height. For single use only.

3.t Metal Ring: Metal rings are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. Metal rings are used in prosthetic restorations together with knob abutment. Metal Rings are designed with 5.0 prosthetic platform diameter and 2.6 mm gingival height. They are for single use only.

3.u Analogs: Analogs are an auxiliary part that enables the patient to restore the chewing function after the Venuscon Dental Implants are surgically placed in the lower or upper jaw bone. Analogs; after the measurements taken with impression pieces are taken on the model; it imitates the implant by mounting on the impression piece. Single use only.

3.v Impression Part: These are the parts used to determine the dimensions of the crown planned to be placed on the implant and superstructure. There are two separate options for classical and digital methods. Single use only.

3.w Parallel Pin: In multiple tooth deficiencies; parallel pins are used to check whether the slots are parallel to each other after the placement of multi unit abutments. In multiple implant placement; the parallelism of the implant sockets is important for the upper prosthesis to be made. It can be used more than once after cleaning and sterilisation.

3.x Depth Gauge: The Depth Gauge is used to check the depth of the created slot after the creation of the implant slots; before the implant is placed in the slot. After measuring the depth, if the neck is suitable, the implant is inserted; if not, the slot is brought to the appropriate size and then the implant is inserted. It can be used more than once after cleaning and sterilisation.

3.y Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool: It is used to ensure that the Screw Retain Drill and Broken Screw Extractor Hand Tool, which are used to remove the broken screw in the mouth, can centre the screw to be removed. The Broken Screw Removal Centring Part is first placed in the area where the screw to be removed and then the screw is removed from the mouth with the Screw Retain Drill or Broken Screw Extractor Hand Tool. It is used for removing broken screws in the patient's mouth. After placement of the Broken Screw Removal Centering Part, the broken screw is removed with the Broken Screw Extractor Hand Tool. It is released as non-sterile. It can be used maximum 15 times after cleaning and sterilisation.

3.z Torque Ratchet / Ratchet : It is used during the placement of implants; during the placement of superstructures with manual parts; for tightening in the mouth according to the recommended torque value. It can be used up to 15 times after cleaning and sterilisation.

3.aa Hex Drivers: It is used for transporting implants or abutments to the patient's mouth and hand tightening with recommended torques. It is recommended not to use more than 15 times. It can be used maximum 15 times after cleaning and sterilisation.

3.bb Hand Tool : During the placement of implants in the mouth; It is used with carrying and clamping parts. It can be used up to 15 times after cleaning and sterilisation.

3.cc Gingival Height Part: The gingival height part is used to measure the level of the gingiva after the implant is placed into the bone. The gingival levels of the gingiva former and abutments to be used are determined after using the gingival height parts. It can be used after cleaning and sterilisation.

4. Product Compatibility

For the compatibility of the products with each other and the products recommended for use, please refer to the additional tables. (You can access the tables from the Medisolaris website. www.medisolaris.com) For implant drill compatibility, please refer to MKK.02 Surgical Instruction for Use.

Annex 1 – MKP.02 – L.03 Medisolaris Implant - Superstructures Compatibility List

Annex 2 – MKP.02 – L.04 Medisolaris Implant - Screws Compatibility List

Annex 3 – MKP.02 – L.05 Medisolaris Implant- Surgical Instruments Compatibility List

Annex 4 – MKP.02 – L.13 Medisolaris Implant-Gingival Height/Impression Part/Scan Body/Screwdriver/Transferring and Tightening Part Compliance List

Annex 5 – MKP.02-L.09 Medisolaris Analogue, Impression Part- Screw Compatibility List

Annex 6 – MKP.02-L.10 Medisolaris Multi Unit Abutment - Scan Body -Impression Part - Burn Out Cap Compatibility List

Annex 7 – MKP.02-L.11 Medisolaris Ratchet-Screwdriver/Transferring and Tightening Part Compatibility List

Annex 8 – MKP.02-L.12 Medisolaris Impression Part-Screw Compatibility List

5. Sterilization and Packaging: Superstructure Components, Healing Heads, Prosthetic Parts and Instruments are packaged unsterilised, sealed in aluminium foil. All system parts must be sterilised before use on patients. All parts must be sterilised by autoclave according to the recommended protocols. The recommended steam sterilisation parameters are as in the table below.

Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacum	Pre-Vacum	3 times < 60 mbar
Temperature	134°C / /270°F	Minimum Drying Time	10 minutes indoor

The packaging of Superstructure Components, Gingiva Formers, Prosthetic Parts, Auxiliary Parts and Hand Tools are not suitable for autoclave sterilisation. The products must be unpacked in a controlled environment before sterilisation, placed in an autoclave bag and then sterilisation must be performed.

6. Placement: Medisolaris dental abutments are placed on the respective implants with carrying and clamping parts and hex keys. Medisolaris ball abutments are self-screwing and are fixed to the implants in the hex slot above the round part using a ball driver. Recommended torque values for screws and superstructures can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants**. Cement should not be used in multi-unit abutment application. Ball abutments and Locator abutments should be used for removable prosthesis for patients who are used to removable prosthesis. Angled options of abutments should be used according to the posture of the implant.

7. Clinical Benefit: Dental implants and Surgical Instruments are a safe, functional alternative treatment. For example, it significantly increases the quality of life of the patient, improves the patient's chewing function, greatly reduces the loss of function, and provides relief from symptoms. It is in the form of natural teeth, provides chewing function, provides a better aesthetic appearance, healthy nutrition and self-confidence with healthy teeth.

A strong and stable dental implant restores the missing tooth so that it looks and functions like a natural tooth. Dental implants are advantageous because the stability of removable prostheses can be improved or the need for fixed tooth-supported or removable prostheses can be met. In a study by Perrotti et al., 2 different implant systems were found to be suitable for long-term treatment, with stable marginal bone levels and a 10-year success rate. Dental implants help you maintain the natural shape of your face and smile. Leaving a gap in your mouth after losing one or more teeth leads to further health problems such as jawbone deterioration and loss. When the natural tooth is not supported, the jawbone deteriorates and loses its strength and stability. Dental implants are the only dental restoration option that protects and stimulates natural bone, helps bone to grow and prevents bone loss.

8. Safety and Performance Information: Test reports that meet the safety and performance criteria of Medisolaris products have been applied and the results of the test reports have been found to be appropriate. Adverse events that may occur in Medisolaris products must be reported by the user or patient to the manufacturer and the relevant competent authority.

Performance Information :

Functioning as a tooth root substitute, Designed to eliminate edentulism in the maxilla and mandible, Providing osseointegration in the jawbone by replacing the tooth/teeth

Safety Information :

Please refer to Contraindication, Residual Risks, Warnings, Possible Side Effects in the instructions for use.

9. Warnings:

The causes of tooth loss should be investigated and the continuity of these causes after implantation should be evaluated. Factors that may affect post-implant healing and long-term use include;

► Oral hygiene habits, ► Plaque control efficiency and cognitive skills specific to the individual, ► Periodontal health of existing teeth,

► Soft tissue health, ► Parafunctions. Venuscon Implants have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment.

Patients with the following diseases are recommended to use implants after being checked by physicians (with consultation if necessary);

- 1) Heart Diseases (using a pacemaker)
- 2) Hypertension
- 3) Diagnosed with or at risk for osteoporosis
- 4) Oral drug use for more than 3 years
- 5) Immune system disorders and immunocompromised patients
- 6) Patients who have received radiotherapy in the past
- 7) Patients with Psychological and Mental Problems
- 8) Smoking - Alcohol and drug users; Patients with teeth clenching and grinding disease.

The location of the implant socket must be correctly adjusted.

There should be at least 1.5 mm between the intact teeth and at least 3 mm between the dental implant. The amount of bone in the vestibule, lingual or palatal should also be at least 1 mm.

The bone temperature should not exceed 47°C during the preparation of the implant bed.

Sterilise surgical instruments before use (according to the manufacturer's recommendation)

Not to use surgical hand tools more than the maximum number of uses

The user is legally capable of performing a surgical operation

Advanced training in the field of pre-operative oral implantation and detailed knowledge of system components

The necessity of giving the Patient Education Brochure to the patient

You should use the prescription prescribed by your dentist for your pain regularly for 2 days

If the pain does not disappear within a few days, the patient should consult a physician.

The sterile sponge placed in the area where dental implant surgery is performed for bleeding control should be bitten with light force for 20 minutes and the sponge should be removed at the end of 20 minutes.

Not eating or drinking anything for two hours after the operation. Do not spit or rinse the mouth with any liquid for two hours.

Avoid eating hard and hot things during the first day (after two hours).

The implanted area should not be touched with the tongue, fingers or a hard object.

Avoid heavy sports and rest after the operation.

Use antiseptic mouthwash for 1 week after the implant application.

Not using tobacco products after implant application

Implants in spoiled sterile packs should not be used

Contaminated implants cannot be sterilised again

Re-use of implants can lead to clinical failure and complications including but not limited to infection, non-integration of implants, bone loss and damage to adjacent structures

Comprehensive patient examination, preoperative imaging, preoperative diagnosis and treatment planning before implant treatment, Inadequate planning and evaluation can lead to implant complications and even implant failure. Immediate and early loading can be achieved with the correct case selection and adherence to the immediate loading protocols established in the literature. If immediate loading is to be performed, the occlusion environment should be checked, multiple implants should be splinted and at least 4 implants in the mandible and 6 implants in the maxilla should be connected to each other with a prosthesis in cases of complete edentulism. Depending on the bone anatomy, the appropriate number, diameter and length of implants should be placed.

Non-axial loads greater than 30° from the long axis of the implant should be avoided. They have not been tested for heating, displacement or image artefacts in the MR environment. The safety of Venuscon in the MR environment is unknown. Scanning a patient with this device may result in patient injury.

10. Information to be given to the patient : MKK.04 Patient Education Brochure must be given to patients by physicians.

11. Disposal: Implants removed from the patient should not be re-implanted as the fatigue status cannot be visually detected. These implants should be considered as biocontaminated and should be treated accordingly. Removed implants should be treated according to their biological hazards and should be disposed of in accordance with the legal regulations of the country and hospital policies.

12. MEANING OF SYMBOLS

	Batch Code		Do not re-use		Sterilized using irradiation
	Catalogue Number		Do not re-sterilize		CE Marking / Notified Body Number
	Used by date		Date of manufacture		Non-Steril
	Do not use if package is damaged and consult instructions for use		Caution		UDI Number
	Consult instructions for use or consult electronic instructions for use		Temperature limit		Manufacturer
	Humidity imitation		Keep away from sunlight		Double sterile barrier System
	Medical Device		Unique Device Identifier		

MEDİSOLARİS SAĞLIK HİZMETLERİ SANAYİ ve TİCARET LİMİTED ŞİRKETİ - Fatih Mahallesi 1191. Sokak No: 28/1 B
Blok Apt. Gaziemir / İZMİR - +90 232 203 6655 / www.medisolaris.com

Ref No: MKK.03

Publish Date: 17.09.2020

Rev. No: 011

Rev Date: 17.10.2024