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MEDISOLARIS Dental Implant Instruction for Use

1-) DENTAL IMPLANTS

Medisolaris (VENUSCON) Dental Implant and Superstructures produced by MEDISOLARIS Sağlık Hizmetleri San. Tic. Ltd. It is and is surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations. VENUSCON Dental Implants are bone level dental implants. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . The target users for Venüscon Dental Implants are dentists.

Superstructures are the device that used to support crowns, bridges or superstructures in the lower and/or upper jaw, used with Dental Implants in completely edentulous or partially edentulous patients.

It is recommended that Medisolaris products not be used for purposes other than their stated intended use.

Medisolaris Dental Implant Basic Information :

For Venuscon Series Implant :

Length : 6 / 8 / 10 / 11.5 / 13 / 16 mm

Diameter : 3.2 / 3.7 / 4.3 / 5.0 mm

Length Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.3	Available	Available	Available	Available	Available	Available
Ø 5.0	Available	Available	Available	Available	Available	Available

1.a Warning : The information presented here is not sufficient for the immediate application of the Medisolaris Dental Implant System. The user must be legally competent to perform surgical operations and must receive advanced training in the field of oral implantation and detailed knowledge of the system components before use.

1.b Legal Liability Disclosure: The use of this product is the sole responsibility of the end user. No liability is assumed for loss or damage caused by misuse of this product, abnormal use (careless handling, sabotage and ignoring safety information), failure to comply with this IFU. This system can only be used with original Medisolaris Dental parts. The use of third party parts renders Medisolaris Sağlık Hizmetleri A.Ş. products ineffective.

1.c Copyright and Patent: The logos of Medisolaris Dental Implants are protected patents and owned exclusively by Medisolaris Sağlık Hizmetleri A.Ş. any use of these marks is prohibited and is subject to prosecution under Turkish and international law.

1.d Education: Dental implant treatment involves complex surgical and restorative procedures. Such procedures should be performed by dentists who have received the necessary training in this field. Inadequate training can lead to surgical and prosthetic complications, resulting in failure of the implant and restoration.

1.e Patient Population : Those with completed bone development (18 years and above), Not pregnant, not in puberty, not breastfeeding, No hard or soft tissue pathology in the implant area, No active periodontal infection or untreated periodontitis, After periodontal treatment, good oral hygiene, good treatment response, non-smoker, being systemically healthy, No unstable endocrine disease (such as diabetes), No local inflammation or mucosal disease (such as lichen planus), Non-osteoporotic, No immunosuppressive disease, Not requiring cancer treatment (chemotherapy/radiotherapy/anti-resorptive drugs), No acquired or inherited coagulation disorder, Not using/has not used IV-based bisphosphonates. No alcohol or drug addiction, No history of allergy to titanium and stainless steel, Not requiring regular medication due to systemic disease, No mental developmental disorder, Patients without recent myocardial infarction or cerebrovascular accident

1.f Raw Materials and Ingredients: Implants and superstructures are made of ISO 5832-3 Ti6AL4V ELI / GRADE23, surgical instruments are made of 1.4116, hand instruments are made of 1.4057 Stainless Steel. The surface of Medisolaris Dental Implant is roughened by sandblasting and passivation SBM method consisting of hydroxyapatite and tricalcium phosphate to obtain surface microstructure, osseointegration, Ra values 2.0 - 2.5 µm.

1.g Indications:

Devices are indicated for total and partially edentulism

Contraindications:

✓ Contraindications ✓ Those whose bone development is not complete (under 18 years of age) ✓ Those whose bone development is not complete (under 18 years of age) ✓ Pregnant, in rapid puberty, or breastfeeding individuals ✓ Presence of hard or soft tissue pathology in the implant area ✓ Individuals with active periodontal infection or untreated periodontitis ✓ Unstable endocrine disease such as uncontrolled diabetes ✓ Presence of local inflammation or mucosal disease such as lichen planus ✓ Individuals diagnosed with osteoporosis	✓ Individuals with immunosuppressive disease ✓ Cancer patients requiring chemotherapy, radiotherapy, or anti-resorptive treatment ✓ Individuals with acquired or hereditary coagulation disorders ✓ Individuals who have used or are using intravenous bisphosphonates ✓ Individuals with alcoholism or substance addiction ✓ Individuals with a history of allergy to titanium or stainless steel ✓ Patients requiring regular medication due to systemic disease ✓ Individuals with mental developmental disorders impairing the ability to maintain oral hygiene ✓ Patients with recent myocardial infarction and cerebrovascular accident
--	---

1.h Residual Risks :

• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening


1.i Potential Side Effects:

• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability

• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening

1.j Placement Configuration of Implants in the Mouth : Since the tooth itself and its root are small in diameter, small diameter implants are placed in areas where the alveolar crest is thinner. After the implant is placed in the alveolar ridge, there should be at least two mm of bone distance around the implant in bucco-lingual and mesio-distal directions. This two mm distance provides the necessary support to the implant during use. In Venuscon series implants; It is recommended that 3.2 mm diameter implants with a narrow platform be placed in the dark area on the side. The image of the oral structure has been prepared in accordance with the tooth and oral region numbering and nomenclature system defined in ISO 3950. The ISO 3950 tooth and oral region numbering and nomenclature system is shown in the image below. In the image on the right, the areas marked in dark color represent the regions in which the use of implants is recommended according to ISO 3950: teeth numbered 12 and 22 in the maxilla, and teeth numbered 42, 41, 31, and 32 in the mandible. The use of implants is not recommended in individuals with incomplete bone development (under 18 years of age) or in cases where deciduous (primary) teeth are present.

It is recommended that the implant size to be used is decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

Likewise, large diameter implants are used in areas where the tooth and bone diameter is thick, that is, there are teeth with more than one root. A large diameter implant is placed in the area where the lower and upper jaw incisors are located, where there is biting force instead of chewing. In Medisolaris implants; Venüscon Ø3.7/4.3/5.0 implants, which have a regular platform, are recommended to be placed in the dark area on the side. According to the ISO 3950 tooth and oral region numbering and nomenclature system provided above, the areas marked in dark color in the image on the side represent the regions in which the use of implants is recommended. These regions correspond to teeth numbered 18, 17, 16, 15, 14, 13, 11, 21, 23, 24, 25, 26, 27, and 28 in the maxilla, and teeth numbered 48, 47, 46, 45, 44, 43, 33, 34, 35, 36, 37, and 38 in the mandible. The use of implants is not recommended in individuals with incomplete bone development (under 18 years of age) or in cases where deciduous (primary) teeth are present. It is recommended that the implant size to be used is decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME	PRODUCT CODE	PRODUCT NAME
SLV 3706	Venuscon Implant Ø3.7mm L:6mm Platform 3.0mm	SLV 4306	Venuscon Implant Ø4.3mm L:6mm Platform 3.0mm
SLV 3708	Venuscon Implant Ø3.7mm L:8mm Platform 3.0mm	SLV 4308	Venuscon Implant Ø4.3mm L:8mm Platform 3.0mm
SLV 3710	Venuscon Implant Ø3.7mm L:11.5mm Platform 3.0mm	SLV 4310	Venuscon Implant Ø4.3mm L:10mm Platform 3.0mm
SLV 3711	Venuscon Implant t Ø3.7mm L:11.5mm Platform 3.0mm	SLV 4311	Venuscon Implant Ø4.3 mm L:11mm Platform 3.0mm
SLV 3713	Venuscon Implant Ø3.7mm L:13mm Platform 3.0 mm	SLV 4313	Venuscon Implant Ø4.3mm L:13mm Platform 3.0mm
SLV 3716	Venuscon Implant Ø3.7mm L:16mm Platform 3.0mm	SLV 4316	Venuscon Implant Ø4.3mm L:1 mm Platform 3.0mm
SLV 5006	Venuscon Implant Ø5.0mm L:6mm Platform 3.0mm	SLV 5011	Venuscon Implant Ø5.0mm L:11.5mm Platform 3.0mm
SLV 5008	Venuscon Implant Ø5.0mm L:8mm Platform 3.0mm	SLV 5013	Venuscon Implant Ø5.0mm L:13mm Platform 3.0mm
SLV 5010	Venuscon Implant Ø5.0mm L:10mm Platform 3.0mm	SLV 5016	Venuscon Implant Ø5.0mm L:16mm Platform 3.0mm

1.k Packaging and Sterilization : Medisolaris Dental implants are packaged sterile with Caps. The products are disposable. Gamma radiation is used to sterilize implants after packaging. Thanks to sealed packaging, implants are prevented from being affected by external environments. For proper storage conditions, the temperature should be between 18-28°C and the humidity should be between 40-60%. Sterilization is valid until the expiry time stated on the label on the package. Implants in damaged sterile packages should not be used. Contaminated implants cannot be resterilized. Reuse of implants may result in clinical failure and complications including, but not limited to, infection, failure to osseointegrate implants, bone loss, and damage to adjacent structures. The implant surface should not be touched by hand, as this may deteriorate its surface properties. All other superstructures and hand instruments should be sterilized by autoclave. For recommended autoclave sterilization parameters, see the Superstructure section and *MKK.02 Surgical Instruction for Use*. Do not use the product whose packaging is damaged.

Implant tube caps are labeled with diameter markings that are color-coded according to the implant diameter.

Ø 3.2 – Turquoise	Ø 3.7 – Green	Ø 4.3 – Purple	Ø 5.0 – Yellow
-------------------	---------------	----------------	----------------

1.l Shelf Life: The expiration date is clearly stated on each implant package.

1.m Treatment Planning: Comprehensive patient examination, pre-operative imaging, pre-operative diagnosis and treatment planning should be done before implant treatment. In patients who cannot achieve the required bone density for the implant, it is recommended to apply for bone augmentation procedures with advanced surgical techniques to increase bone density before the application and then perform the application. Inadequate planning and evaluation can lead to implant complications or even implant failure.

1.n Preparation of the Surgical Area: Preparation of the implant area should be done with the surgical drills in the Medisolaris Dental surgical set. The drilling speed should be decided by the physician according to the hardness level of the bone. Drilling should be performed under abundant irrigation. Care should be taken to prevent the bone from getting too hot. Bone temperature should not exceed the upper limit of 47 °C. The use of surgical drills that have lost their sharpness should be avoided. Because these may cause unnecessary trauma to the bone. It is recommended to follow the drilling sequence according to the implant diameter to be used. For detailed information about the surgical operation, see *MKK.02 Surgical Instruction for Use*.

1.o Surgical Placement of Implants: Check the length and diameter of the implant on the product label. Attach the (3) labels included in the packaging to the patient's file. Open the outer cardboard package in a controlled environment and remove the thermoformed packaging. Open the tyvek paper on the thermoform and empty the sterile tube into the surgical area. The implant tube consists of a single part: Open the implant tube by turning the top cap. Use the implant carrier in the set to move the implant into the surgically prepared implant slot. The ratchets or physiodispenser included in the set should be used to place the implant. Medisolaris Covers should be closed with the help of a 1.25 mm hex key. Caps and healing caps should be torqued by hand.

Recommended Torq Values

Process	Recommended Torq Values	Process	Recommended Torq Values
Placement of Implants	35-40 Ncm	Screw Venuscon	30 Ncm
Gingiva Former	15 Ncm	Multi-Unit occlusal	20 Ncm
Cover	15 Ncm	Self-Screwing Abutments	30 Ncm

For detailed information see *MKK.02 Surgical Instruction for Use*. *Recommended used at 700 RPM.

1.p Healing Process of Medisolaris Dental Implants: In general, implants should be left to heal for 2 to 4 months, depending on the patient's bone quality, type and general health status. Immediate loading can be performed with the right case selection and adherence to the immediate loading protocols accepted in the literature. If immediate loading is to be performed, the occlusion environment should be checked. If there are multiple implants, they should be splinted and in cases of complete edentulism, at least 4 implants in the mandible and 6 implants in the maxilla should be connected to each other with a prosthesis. Depending on the bone anatomy, the appropriate number, diameter and length of implants should be placed.

1.q Information to be communicated to the patient after implant placement

It is important that patients follow the instructions given below after dental implant surgery to protect their implants for life. For detailed information, please refer to MKK.04 Patient Education Brochure. It is recommended to be given to patients after the operation.

In dental implant surgical treatment, it is recommended to wait at least 4 months in the upper jaw and at least 3 months in the lower jaw for the implant to fuse with the jaw bone.

After the dental implant is placed, the patient will feel some pain and discomfort. The pain level is usually mild. After the operation, the dentist may prescribe painkillers. With painkillers, the pain can be controlled. These should be taken regularly for 2 days. The first dose of painkillers must be taken before the local anaesthetic wears off. The pain will begin to subside within a few days. If the pain does not go away after a few days, you should consult your doctor.

1. After the dental implant, the dentist may prescribe antibiotics and antiseptic mouthwash to the patient if deemed appropriate.
2. A sterile sponge is placed on the area where dental implant surgery is performed to control bleeding. The sterile sponge should be bitten by the patient with light force for 20 minutes. At the end of 20 minutes, the sponge should be removed.
3. The patient should not eat or drink anything for two hours after the operation. The patient should not spit or rinse his/her mouth with any liquid for two hours. After two hours, the patient should not eat hard and hot foods on the first day. Hot foods increase both bleeding and the risk of infection.
4. It is normal to have bleeding in the form of leakage during the first 24 hours after the procedure. To control bleeding, the patient bites gauze for 30 minutes. If the bleeding continues for more than 2 days, you should consult our physician. Aspirin and aspirin compounds should not be taken as painkillers because they cause bleeding tendency.
5. Swelling may occur after implant surgery. This is a normal condition. If the pain and swelling do not disappear after 5-6 days, you should consult your physician.
6. The patient should avoid tampering with the area where the dental implant is placed with the tongue, finger or hard instruments.
7. The patient should rest a little after the implant surgery and avoid strenuous exercise for 48 hours.
8. If the patient is using removable dentures, he/she should not use them until the date that the dentist deems appropriate.
9. In dental implant surgery, the wound surfaces are closed with sutures. The sutures must be removed by the dentist after one week at the earliest. The patient is informed by the dentist about the early period in which the treated tooth part can be used.
10. After dental implant surgery, the patient should continue oral care after the first 24 hours. The patient should use antiseptic mouthwash prescribed by the dentist three times a day for the first week after the surgery. Antiseptic mouthwash should not be used after one week. After antiseptic mouthwash, the patient can gargle with warm salt water. Add 1 teaspoon of salt to 1 glass of water and gargle 3-4 times a day to disinfect the area. This will help clean the healing area.
11. After 48 hours, the surgical site should be lightly brushed. Other teeth can be brushed and cleaned normally after the operation. There should be no food residue in the surgical site and the patient should rinse the mouth with warm water after meals.
12. Smoking delays the healing process and prevents dental implants from staying in the mouth for a long time, so the patient should stop smoking cigarettes and similar tobacco before starting dental implant treatment.
13. It is very important to follow all these instructions mentioned above for the success of dental implant surgery. If there are signs of infection after implant surgery or if the patient experiences nausea, fever or vomiting, he/she should consult his/her dentist.

1.r Long-term care of dental implants: In order to be able to use implants for a lifetime without any problems, it is important that they are properly cared for as described above.

1.s The importance of oral hygiene in dental implant patients: Implants are more susceptible to plaque and calculus build-up, inflammation and bone loss than natural teeth. In implant patients, good oral hygiene is very important. The patient should brush teeth correctly and clean the interdental and interimplant areas correctly. Oral hygiene kits are available for the private use of implant patients, which can be used by the patient. If the patient follows the oral hygiene instructions, inflammation of the periimplant tissues will be minimal. Removal of plaque and calculus deposits from the teeth will reduce the risk of periimplantitis, the most common cause of implant failure.

1.t Regular dental visits and examinations for implant patients: In patients treated with dental implant prostheses, check-ups are usually recommended every 3 months for the first 1-2 years and then every 6 months thereafter. At each visit, oral hygiene instructions are reinforced and the patient's ability to maintain good oral hygiene is checked. Patients with implant-retained removable prostheses are given some special instructions. The patient is told how to insert and remove the prosthesis without causing any damage to the implant or prosthesis. Excessive pressure should not be applied to the restoration placed on the implant. Restorations on implants are not as strong as natural teeth. The patient should avoid eating hard foods that can break the restoration. Avoid irritation of the implant area.

1.u Lifespan of dental implants: The lifespan of dental implants depends on many factors, such as patient health and the care provided by the patient. The lifespan of implants depends on how old the patient had the implants, how they live and how they care for them. For the success of the implants, the patient should maintain good oral hygiene and visit the dentist regularly for check-ups. Although a clear lifetime cannot be defined for the device, it is thought to have a lifetime of 10 years under normal conditions. Implants have a shelf life of 5 years and superstructures do not have a shelf life.

Implants are an effective solution for missing teeth and can be used for a lifetime with good care. Brushing, flossing and regular visits to the dentist after the implant prosthesis are important to prolong the implant's stay in the mouth. Implants can replace a single tooth, a group of teeth or all upper and lower jaw teeth.

1.v Opening Sterile Packages: To take Medisolaris Dental implants from the product packaging with an angldurva, implant keys suitable for angldurva are used. To take the implants manually, the implant key suitable for manual removal is used. When taking the implants, the size suitable for the area to be placed in the mouth should be selected. For all diameters of the Venuscon implant, the 2.3 mm hexagon crimping spanner included in the surgical set should be used.

In order to attach the closing screws of Medisolaris Dental implants to the implants, the caps that are compatible in diameter that come out of the product are placed with a hex key.

STEP 1: Open the implant box and remove the sterile packaged thermoform.

STEP 2: The back of the thermoform is opened. If the adhesion area is torn, broken or opened, the product should not be used. The implant should be opened in a sterile environment where the surgical operation will be performed.

STEP 3: If the implant is to be taken with an angldurva; the key suitable for the angldurva is attached to the tip of the angldurva. If the implant is to be taken manually, the manual implant key is attached to the torque racket.

STEP 4: The cap of the tube is held upwards. After the cap of the tube is opened by turning the cap, the implant in the sling is removed with an angldurva implant key or manual implant key.

STEP 5: To remove the implant from the sling, the angldurva implant key or manual implant key is inserted into the implant. Make sure that the keys are fully inserted into the implant so that the implant does not fall out.

STEP 6: When the implants are removed from the discs in the sling, they are removed by moving back and forth.

STEP 7: If an angldurva is used to place the implant, the appropriate torque is adjusted from the motor and placed in the slot opened in the bone. If the implant is placed manually, the implant is removed from the slot and placed in the opened slot after the appropriate torque is adjusted before taking the implant to the torque raschette.

STEP 8: In order to place a cover on the implant, the cover in the tube is removed with the help of a key suitable for the cover in the surgical operation room.

STEP 9: The cap removed from the tube is placed on the implant without manual contact.

If the products are not to be taken directly with hand instruments, they should be opened on the sterile cloth.

1.y Post-Operative Care: Close follow-up of the patient after surgery is very important. It is recommended to use appropriate oral rinses and to ensure good oral care and not to use tobacco products.

2-) Medisolaris Dental Cover

The closure screw is used to cover and protect the implant after placing the implant in the bone. The closure screw can remain in the mouth for 2-3 months during the surgical healing period.

2.a Raw Material : ISO 5832-23 Ti6AL4V ELI

2.b Indications: The caps come out of the implant package used. It is used in two-stage implant surgery to prevent soft or hard tissue growth on the hollow inner surface of the implant.

2.c Contraindications: Contraindications are mentioned above.

2.d Possible Side Effects: Possible side effects are mentioned above.

2.e Placement: The caps are torqued by hand using a 1.25 mm hex spanner. Recommended torque values can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants**.

2.f Packaging and Sterilisation: The caps are packaged together with the implant in a sterile tube. If the package is intact, sterility is valid until the expiration date indicated on the package. In cases where the cap is packaged individually, the product is single-packed without being sterile and wrapped in aluminum foil; it must be sterilized before use by the autoclave method according to the parameters specified under Section 5, **Sterilization and Packaging**, of the Instructions for Use.

The packaging of the caps is not suitable for autoclave sterilization. Prior to sterilization, the products must be removed from their packaging in a controlled environment, placed into packaging suitable for sterilization, and then sterilized. The calibration of the autoclave must be performed at least once a year, and the results must be within acceptable limits.

3-) Medisolaris Dental Superstructures, Gingiva Formers, Prosthetic Parts and Instruments

Superstructures are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with dental implants in completely edentulous or partially edentulous patients.

3.a Raw Material : ISO 5832-3 / Ti6AL4V ELI (Grade 23) - Stainless Steel 1.4116 and 1.4057

3.b Gingiva Formers: Healing caps are used after the healing of the implant in the jawbone is completed and the overlying soft tissue is removed, until prosthetic treatment for gum shaping. Healing caps remain in the mouth for 15 days or less. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are for single use only.

3.c Temporary Abutment: Temporary abutment is used with Medisolaris implant in completely edentulous or partially edentulous patients. In immediate loading treatments, they are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw for more than 24 hours. Temporary abutment is an endosseous dental implant abutment designed for single use as a temporary prosthesis during the healing process while the permanent prosthesis is being produced. Temporary abutments are compatible with all diameters of Medisolaris implant for Venuscon series. It is fixed to the implant with the help of a screw. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are for single use only.

3.d Straight Abutment: In order to use a cementable abutment, the implant must be placed in the patient without an angle. Cementable Abutments are available for both brands. It is compatible with all diameters of Medisolaris Dental implants for Venuscon series. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is bonded on it with the help of cement. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.e Angled Abutment: Angled abutments facilitate prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid anatomical structure when placing the implant. The angled abutment is compatible with all diameters of the Medisolaris implant for the Venuscon series. It is fixed to the implant with the help of a screw. The angled abutment is used in cemented restorations. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.f Anatomical Straight Abutment: The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This offers a more aesthetic option to the physician. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it emerges through the gum. In a possible gum recession, the underlying titanium is prevented from being seen. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zirconia is glued on it. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is disposable.

3.g Anatomically Angled Abutment : The use of an anatomically angled abutment facilitates prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid the anatomical structure when placing the implant. In addition, the use of aesthetic angled abutment shortens the healing time. The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This gives the dentist a more aesthetic option. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it comes out of the gum. In a possible gum recession, the underlying titanium is prevented from being seen. For each brand (Medisolaris Dental) implant, the same brand Anatomically Angled Abutment should be used. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is glued on it. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.h Multi Unit Abument: It is the external connection of the implant embedded flat into the bone. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different gum options are available. It is self-screwed and fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is disposable.

3.i Multi Unit Angled Abutment: It is the external connection of the implant embedded into the bone at an angle. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different angle and gum options are available. It is fixed to the implant with a screw. A tooth is made of materials such as porcelain or zircon. There are two types of multi-unit angled abutment compatible with Medisolaris implants: hex and non-hex. The hex structure is recommended to be used in single tooth restorations, while the non-hex structure is recommended to be used in bridge, bar and hybrid treatments. Hex structure is designed to prevent rotation in single tooth applications. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.j Ball Abutment: It is the external connection of the implant embedded flat into the bone. It is used in cases where the physician wants to apply removable prosthesis to the patient. Different diameters and gum lengths are available. It is self-screwed and fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.k Tibase Cad-Cam Abutment: It is the platform of the abutment inside the implant. Zirconia or custom (personalised) abutment is produced on it. It is produced for the purpose of using a company-specific platform. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.l Locater Abutment: It is an abutment system with low vertical height used for implant over implant removable prosthesis in completely edentulous jaws. The locater abutment piece is fixed to the implant with a metal cap and rubber in the removable prosthesis. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.m Universal Straight Abutment: Universal Straight Abutment are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. It is milled and made suitable for the patient. It is fixed to the implant with the help of a screw. When preparing the universal abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant completely or the connection will not be provided. After the universal abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the universal abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will

decrease. The milled universal abutments should be washed at 55 - 60 °C for 15 minutes. The products should not be sterilised and used without washing. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.n Premill Abutment: It is produced for CAD-CAM devices. The laboratory connects this superstructure to the device and makes a patient-specific abutment. While preparing the premill abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant properly or the connection will not be provided. After the premill abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the premill abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. Premilled abutments that have undergone a milling process must be washed at 55–60 °C for 15 minutes. The products must not be sterilized or used without first performing the washing process. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are single-use only.

3.o Ti-base Casting Abutment: For use in personalised cast prosthetic restorations in single or multiple implants. It is fixed to the implant with the help of a screw. Implants for Venuscon series in Medisolaris Dental implants are designed to be compatible with all superstructures. Interface connection type is the same. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.p Screws: The screw allows the superstructures to be mounted to the implant. The screws used with Measuring Parts and Analogues are different. It is packaged and sold together with the related product. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.q Titanium Cap: In completely or partially edentulous patients; These are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants. Titanium Caps are designed as Hex and No-Hex. Both caps are compatible with Venuscon Multi Unit Abutments. Ø 5.0 prosthetic platform diameter and 12.0 mm gingival height. No coating process is applied on the surface. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.r Multi-Unit Tibase Abutment: Multi-Unit Tibase abutment; These are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants in completely edentulous patients. Multi-Unit Tibase abutment is designed as hex and no-hex. Both abutments are compatible with Venuscon Implants. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.s Comfort Cap: The comfort cap prevents the patient from being uncomfortable by covering the multi unit abutment in the patient's mouth until the crown is loaded in Multi unit abutments. Comfort Caps are designed with 5.0 prosthetic platform diameter and 4.8 mm gingival height. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . For single use only.

3.t Metal Ring: Metal rings are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. Metal rings are used in prosthetic restorations together with knob abutment. Metal Rings are designed with 5.0 prosthetic platform diameter and 2.6 mm gingival height. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are for single use only.

3.u Analogs: Analogs are an auxiliary part that enables the patient to restore the chewing function after the Venuscon Dental Implants are surgically placed in the lower or upper jaw bone. Analogs; after the measurements taken with impression pieces are taken on the model; it imitates the implant by mounting on the impression piece. Manufactured from stainless steel 1.4057 material .Single use only.

3.v Impression Parts: These are the parts used to determine the dimensions of the crown planned to be placed on the implant and superstructure. There are two separate options for classical and digital methods. Manufactured from Ti6AL4V ELI and stainless steel 1.4057 material . Single use only.

3.w Parallel Pin: In multiple tooth deficiencies; parallel pins are used to check whether the slots are parallel to each other after the placement of multi unit abutments. In multiple implant placement; the parallelism of the implant sockets is important for the upper prosthesis to be made. It can be used more than once after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material

3.x Depth Gauge: The Depth Gauge is used to check the depth of the created slot after the creation of the implant slots; before the implant is placed in the slot. After measuring the depth, if the neck is suitable, the implant is inserted; if not, the slot is brought to the appropriate size and then the implant is inserted. It can be used more than once after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material

3.y Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool:

It is used to ensure that the Screw Retain Drill and Broken Screw Extractor Hand Tool, which are used to remove the broken screw in the mouth, can centre the screw to be removed. The Broken Screw Removal Centring Part is first placed in the area where the screw to be removed and then the screw is removed from the mouth with the Screw Retain Drill or Broken Screw Extractor Hand Tool. It is used for removing broken screws in the patient's mouth. After placement of the Broken Screw Removal Centering Part, the broken screw is removed with the Broken Screw Extractor Hand Tool. It is released as non-sterile. It can be used maximum 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material

3.z Torque Ratchet / Ratchet : It is used during the placement of implants; during the placement of superstructures with manual parts; for tightening in the mouth according to the recommended torque value. It can be used up to 15 times after cleaning and sterilisation. Manufactured from Ti6AL4V ELI (Grade 23) and Stainless Steel 1.4057 material .

3.aa Hex Drivers: It is used for transporting implants or abutments to the patient's mouth and hand tightening with recommended torques. It is recommended not to use more than 15 times. It can be used maximum 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4116 and stainless steel 1.4057 .

3.ab Hand Tools : During the placement of implants in the mouth; It is used with carrying and clamping parts. It can be used up to 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4116 and stainless steel 1.4057 .

4. Product Compatibility

For the compatibility of the products with each other and the products recommended for use, please refer to the additional tables. (You can access the tables from the Medisolaris website. www.medisolaris.com) For implant drill compatibility, please refer to MKK.02 Surgical Instruction for Use.

Annex 1 – MKP.02 – L.03 Medisolaris Implant - Superstructures Compatibility List

Annex 2 – MKP.02 – L.04 Medisolaris Implant - Screws Compatibility List

Annex 3 – MKP.02 – L.05 Medisolaris Implant- Surgical Instruments Compatibility List

Annex 4 – MKP.02 – L.13 Medisolaris Implant-Gingival Height/Impression Part/Scan Body/Screwdriver/Transferring and Tightening Part Compliance List

Annex 5 – MKP.02-L.09 Medisolaris Analogue, Impression Part- Screw Compatibility List

Annex 6 – MKP.02-L.10 Medisolaris Multi Unit Abutment - Scan Body -Impression Part - Burn Out Cap Compatibility List

Annex 7 – MKP.02-L.11 Medisolaris Ratchet-Screwdriver/Transferring and Tightening Part Compatibility List

Annex 8 – MKP.02-L.12 Medisolaris Impression Part-Screw Compatibility List

5. Sterilization and Packaging: Superstructure Components, Gingiva Formers, Prosthetic Parts and Instruments are packaged unsterilised, sealed in aluminium foil. All system parts must be sterilised before use on patients. All parts must be sterilised by autoclave according to the recommended protocols. The recommended steam sterilisation parameters are as in the table below.

Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacum	Pre-Vacum	3 times < 60 mbar
Temperature	134°C / /270°F	Minimum Drying Time	10 minutes indoor

The packaging of Superstructure Components, Gingiva Formers, Prosthetic Parts, Auxiliary Parts and Hand Instruments are not suitable for autoclave sterilisation. The products must be unpacked in a controlled environment before sterilisation, placed in an autoclave bag and then sterilisation must be performed. The calibration of the autoclave must be performed at least once a year, and the results must be within acceptable limits.

6. Placement: Medisolaris dental abutments are placed on the respective implants with carrying and clamping parts and hex keys. Medisolaris ball abutments are self-screwing and are fixed to the implants in the hex slot above the round part using a ball driver. Recommended torque values for screws and superstructures can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants.**

Cement should not be used in multi-unit abutment application. Ball abutments and Locator abutments should be used for removable prosthesis for patients who are used to removable prosthesis. Angled options of abutments should be used according to the posture of the implant.

7. Clinical Benefit: Dental implants and Surgical Instruments are a safe, functional alternative treatment. For example, it significantly increases the quality of life of the patient, improves the patient's chewing function, greatly reduces the loss of function, and provides relief from symptoms. It is in the form of natural teeth, provides chewing function, provides a better aesthetic appearance, healthy nutrition and self-confidence with healthy teeth.

A strong and stable dental implant restores the missing tooth so that it looks and functions like a natural tooth. Dental implants are advantageous because the stability of removable prostheses can be improved or the need for fixed tooth-supported or removable prostheses can be met. In a study by Perrotti et al., 2 different implant systems were found to be suitable for long-term treatment, with stable marginal bone levels and a 10-year success rate. Dental implants help you maintain the natural shape of your face and smile. Leaving a gap in your mouth after losing one or more teeth leads to further health problems such as jawbone deterioration and loss. When the natural tooth is not supported, the jawbone deteriorates and loses its strength and stability. Dental implants are the only dental restoration option that protects and stimulates natural bone, helps bone to grow and prevents bone loss.

8. Safety and Performance Information: Test reports that meet the safety and performance criteria of Medisolaris products have been applied and the results of the test reports have been found to be appropriate. Adverse events that may occur in Medisolaris products must be reported by the user or patient to the manufacturer and the relevant competent authority.

Performance Information :

Functioning as a tooth root substitute, Designed to eliminate edentulism in the maxilla and mandible, Providing osseointegration in the jawbone by replacing the tooth/teeth

Safety Information :

Please refer to Contraindication, Residual Risks, Warnings, Possible Side Effects in the instructions for use.

9. Warnings: The causes of tooth loss should be investigated and the continuity of these causes after implantation should be evaluated. Factors that may affect post-implant healing and long-term use include;

► Oral hygiene habits, ► Plaque control efficiency and cognitive skills specific to the individual, ► Periodontal health of existing teeth,

► Soft tissue health, ► Parafunctions.

Patients with the following diseases are recommended to use implants after being checked by physicians (with consultation if necessary);

1. In patients with heart disease (those using a pacemaker), implant treatment is recommended only after a cardiology consultation. (Subject to the physician's evaluation.)
2. In patients with hypertension, blood pressure should be measured and evaluated before implant treatment. The patient should be informed about the regular use of antihypertensive medications. (Subject to the physician's additional evaluation.)
3. Patients who have not been diagnosed with osteoporosis but are at risk: implant treatment should be performed after an osteoporosis consultation. (Subject to the physician's additional evaluation.)
4. In patients who have used oral medications for more than 3 years, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
5. In patients with immune system disorders or immunosuppressed patients, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
6. In patients who have previously received radiotherapy or chemotherapy, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
7. In patients with psychological and mental disorders, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
8. Patients who smoke, consume alcohol, or use medications; and patients with bruxism (teeth clenching or grinding).

The implant site must be positioned correctly.

The distance from adjacent natural teeth must be at least 1.5 mm, and the distance from another dental implant must be at least 3 mm.

The amount of bone on the vestibular, lingual, or palatal side must be at least 1 mm.

During preparation of the implant bed, the bone temperature must not exceed 47 °C.

Surgical instruments must be sterilized before use (according to the manufacturer's recommendations).

Surgical instruments must not be used beyond their maximum number of uses.

The user must be legally qualified to perform surgical procedures.

Before use, detailed information about oral implantology and system components must be obtained through advanced training.

The patient must be provided with a Patient Information Brochure.

For pain management, the prescription written by the dentist should be used regularly for the duration specified by the physician.

If the patient's pain does not subside within a few days, they must consult the physician.

The sterile gauze placed at the surgical site for bleeding control after dental implant surgery should be bitten with light pressure for 20 minutes, unless otherwise instructed by the physician, and removed after 20 minutes.

For two hours after surgery, the patient should not eat or drink anything. During these two hours (unless otherwise instructed by the physician), the patient should not spit or rinse the mouth with any liquid. For the rest of the first day (after the initial two hours), hard and hot foods should be avoided.

The implant site must not be touched with the tongue, fingers, or any hard object.

After the operation, heavy sports should be avoided and the patient should rest.

An antiseptic mouthwash should be used for one week after implant treatment.

Tobacco products must not be used after implant treatment.

Implants in damaged sterile packages must not be used.

Contaminated implants must not be re-sterilized.

Reusing implants may lead to clinical failure and complications including, but not limited to, infection, lack of implant integration, bone loss, and damage to adjacent structures.

A comprehensive patient examination, preoperative imaging, preoperative diagnosis, and treatment planning must be performed before implant treatment. Inadequate planning and evaluation may lead to implant complications or even implant failure.

MR-Conditional

Non-clinical testing has demonstrated that implants belonging to the Medisolaris Venuscon and Merkurhex product families, when used in combination with straight and angled abutments, can be conditionally scanned in whole-body MRI systems with a 3-Tesla static magnetic field. A patient with these devices may undergo an MRI scan under the following conditions:

- **Static Magnetic Field:** 3-Tesla
- **Maximum SAR Value:** Whole-body averaged Specific Absorption Rate (SAR) reported by the scanner of 2 W/kg
- **Magnetically Induced Displacement Angle:** 2°
- **Scan Duration:** 15 minutes of continuous scanning

- **RF-Induced Heating:** Under the scanning conditions specified above, the maximum temperature increase generated by the devices in the worst-case scenario is 2.6°C for the Venuscon family and 2.4°C for the Merkurhex family
- **Image Artifact:** The quality of the MR image may be affected by the presence of the implant. In testing, the maximum artifact extent was measured as 11 mm for Spin Echo sequences and 19 mm for Gradient Echo sequences. Therefore, if the imaging region is close to the implant location, optimization of MRI parameters is recommended
- **Mechanical Safety:** It has been verified that the magnetically induced displacement force and torque are lower than the gravitational forces encountered during normal daily activities and do not pose a clinical risk

Caution: This safety information applies only to the tested Medisolaris implant and abutment combinations; MR safety cannot be guaranteed if the device is used in combination with components from other manufacturers.

10. Information to be given to the patient : MKK.04 Patient Education Brochure must be given to patients by physicians.

11. Disposal: After the implant removal procedure, the implant component may be infectious, and the instruments may be sharp and piercing, so they are classified as “medical waste.” Implants are classified as “infectious waste”; instruments used in implant placement (e.g., drills, clamping devices) are classified as “sharp waste” due to their cutting and piercing properties.

These containers are manufactured from medium-density polyethylene raw material that is resistant to tearing, puncturing, bursting, and transport, and are leak-proof, double-bottom stitched, and gusset-free, with a double-layer thickness of 100 microns, a minimum lifting capacity of 10 kilograms, and bearing the black “International Biohazard” emblem and the red “CAUTION! MEDICAL WASTE” label in a visible size on both sides. The bags are filled to no more than ¾ capacity, tightly sealed, and, where necessary, each bag is placed inside another bag with the same characteristics to ensure complete leak-proofing. These bags cannot be recycled or reused in any way. The contents of medical waste bags must not be compressed in any way; medical waste must not be removed from the bag, emptied, or transferred to another container.

Sharps and piercing medical waste must be collected separately from other medical waste in containers or boxes made of plastic or laminated cardboard with the same properties, bearing the black “International Biohazard” symbol and the words “CAUTION! CUTTING AND PUNCTURING MEDICAL WASTE” in black letters. These collection containers shall be filled to no more than ¾ capacity, sealed, and placed in medical waste bags. Once sharps containers are full, they must not be compressed, opened, emptied, or recycled. Implant waste must not be mixed with household waste under any circumstances.

Medical waste collected by hospitals, clinics, dental centers, relevant municipal departments, or companies authorized by municipalities is disposed of by the relevant municipality or authorized medical waste disposal facilities.

Personal protective equipment (gloves, masks, goggles, aprons) must be used before and during disposal. Disposal operations must be carried out in accordance with the Medical Waste Regulation.

Each country must act in accordance with its waste management legislation.

12. SSCP

Following the certification of Medisolaris products under MDR 2017/745, the SSCP link will be added.

13. MEANING OF SYMBOLS

	Batch Code (All Products)		Do not re-use (Dental Implant – Abutment – Gingiva Former – Covers)		Sterilized using irradiation (Dental Implants)
	Catalogue Number (All Products)		Do not re-sterilize (Dental Implant – Abutment – Gingiva Former – Covers)		CE Marking / Notified Body number (All Products)
	Used by date (All Products)		Date of manufacture (All Products)		Non-Steril (All product expect dental implants)
	Do not use if package is damaged and consult instructions for use (All Products)		Caution (All Products)	 (01)01234567891011 (11)200622 (17)240622 (21)88888888	UDI Number (All Products)
	Consult electronic instructions for use (All Products)		Temperature limit (All Products)		Manufacturer (All Products)

	Humidity imitation (All Products)		Keep away from sunlight (All Products)		Double sterile barrier System (Dental Implant)
	Medical Device (All Products)		Unique Device Identifier (All Products)		Keep Dry (All Products)
	MRI Conditional (Dental Implant – Abutment – Screw)				

MEDİSOLARİS SAĞLIK HİZMETLERİ SANAYİ ve TİCARET LİMİTED ŞİRKETİ - Fatih Mahallesi 1191. Sokak No: 28/1 B Blok
Apt. Gaziemir / İZMİR - +90 232 203 6655 / www.medisolaris.com

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www.medisolaris.com

A. Medisolaris Dental Implants

There are two separate series of Medisolaris Dental Implants, Merkurhex and Venuscon. Medisolaris Dental Implants are bone-level dental implants that are surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations.

A.1 Medisolaris Dental Implants Features

Medisolaris Dental Implant System is used for the treatment of mandibular and maxillary; single and partial edentulism.

A device specially designed to be surgically implanted in or over the mandibular or maxillary bone as a means of support and retention against dislodgement of the dental prosthesis. It is either threaded through the gum by the dentist (part of the implant emerges through the gum for direct support) or completely embedded under the gum (for removable denture support only).

Medisolaris Dental Implant is a bone-level dental implant that is surgically placed in the upper or lower jaw bone to support artificial dentures used to restore the patient's chewing function. From top to bottom, the implant's core geometry descends conically. Implants with a tapered core have higher tensile force and insertion torque. Reducing the core diameter without increasing the outer diameter too much increases the contact area between the bone implant and increases the tensile strength.

Diameter and Length Options

- Merkurhex Implant

Lenght Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.2	Available	Available	Available	Available	Available	Available
Ø 4.7	Available	Available	Available	Available	Available	Available
Ø 5.2	Available	Available	Available	Available	Available	Available

- Venuscon Implants

Lenght Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.3	Available	Available	Available	Available	Available	Available
Ø 5.0	Available	Available	Available	Available	Available	Available

1.a Warnings : The information presented in the Instruction for Use is not sufficient for the immediate application of the Medisolaris Dental Implant System. The user must be legally qualified to perform surgical operations and must have advanced training in oral implantation and detailed knowledge of the system components before use.

Surgical instruments must be sterilized (according to the manufacturer's recommendation) before use.

Surgical instruments should not be used more than the maximum number of times, and the location of the implant socket should be adjusted correctly.

There should be at least 1.5 mm between the intact teeth and at least 3 mm with the dental implant. The amount of bone in the vestibule, lingual or palatal should also be at least 1 mm.

The bone temperature should not exceed 47C° during the preparation of the implant bed.

1.b Legal Liability Disclosure: The use of this product is the sole responsibility of the end user. No liability is assumed for loss or damage caused by misuse of this product, abnormal use (careless handling, sabotage and ignoring safety information), failure to comply with this IFU. This system can only be used with original Medisolaris Dental parts. The use of third party parts renders Medisolaris Sağlık Hizmetleri A.Ş. products ineffective.

1.c Copyright and Patent: The logos of Medisolaris Dental Implants are protected patents and owned exclusively by Medisolaris Sağlık Hizmetleri A.Ş. any use of these marks is prohibited and is subject to prosecution under Turkish and international law.

1.d Education: Dental implant treatment involves complex surgical and restorative procedures. Such procedures should be performed by dentists who have received the necessary training in this field. Inadequate training can lead to surgical and prosthetic complications, resulting in failure of the implant and restoration.

1.e Patient Populations:

- Those with completed bone development (18 years and above)
- Not pregnant, not in puberty, not breastfeeding
- No hard or soft tissue pathology in the implant area
- No active periodontal infection or untreated periodontitis
- After periodontal treatment: good oral hygiene, good treatment response, non-smoker, systemically healthy
- No unstable endocrine disease (e.g., diabetes)
- No local inflammation or mucosal disease (e.g., lichen planus)
- Non-osteoporotic
- No immunosuppressive disease
- Not requiring cancer treatment (chemotherapy/radiotherapy/anti-resorptive drugs)
- No acquired or inherited coagulation disorder
- Not using / has not used IV-based bisphosphonates
- No alcohol or drug addiction
- No history of allergy to titanium or stainless steel
- Not requiring regular medication due to systemic disease
- No mental developmental disorder
- Patients without recent myocardial infarction or cerebrovascular accident

1.f Raw Material : Stainless Steel and titanium.

1.g Contraindications:

- Those whose bone development is not complete (under 18 years of age)

- Pregnant, in puberty, or breastfeeding individuals
- Individuals with hard or soft tissue pathology in the implant area
- Individuals with active periodontal infection or untreated periodontitis
- Individuals with poor oral hygiene after periodontal treatment, inadequate treatment response, smoking habits, or systemic disease
- Individuals with unstable endocrine disease (such as uncontrolled diabetes)
- Individuals with local inflammation or mucosal disease such as lichen planus
- Individuals diagnosed with osteoporosis
- Individuals with immunosuppressive disease
- Cancer patients requiring chemotherapy, radiotherapy, or anti-resorptive treatment
- Individuals with acquired or hereditary coagulation disorders
- Individuals who have used or are using intravenous bisphosphonates
- Individuals with alcohol or substance addiction
- Individuals with a history of allergy to titanium or stainless steel
- Patients requiring regular medication due to systemic disease
- Individuals with mental developmental disorders impairing the ability to maintain oral hygiene
- Patients with recent myocardial infarction or cerebrovascular accident



1.h Residual Risks : Infection, injury during surgical procedure, perforation and parasthesia after surgical intervention, inability to regain chewing function or delay in regaining chewing function, cutting the gum and hurting the patient, swallowing the implant, swallowing the Medisolaris drill extensions, swallowing the Medisolaris drill, swallowing the clamping piece, the need for a re-operation, Failure to create an implant socket, resulting in the product not serving the purpose of an artificial tooth root, creating stress in the bone - necrosis (dead) bone formation and failure to ensure osseointegration of the implant, prolongation of the operation, Periimplantitis (Implant), Perimucositis (Implant), Implant mobility (Implant), Failure to achieve osseointegration (Implant), Primary Stability (Implant), Inappropriate level of stable bone and soft tissue (Implant)

1.i Potential Side Effects:

• Infection	• Difficulty opening the mouth
• Bleeding and hematoma	• Tissue numbness
• Injuries to the nerves, sinus and nasal cavity	• Tingling in tissues
• Incorrect positioning of the implant	• Sensory Alteration
• Incision line patency	• Damage to neighboring tooth
• Bruising	• Sinus and nasal floor injury in the upper jaw
• Swelling	

1.j Placement Configuration of Implants in the Mouth:

Small diameter implants are placed in areas where the alveolar crest is thinner because the tooth itself and the root are small in diameter. After placing an implant in the alveolar crest, at least two mm of bone distance should remain around the implant in the bucco-lingual and medio-distal directions. This two mm distance provides the necessary support to the implant during the use of the implant. In Venuscon and Merkurhex series implants; 3.2 mm diameter implants with a narrow platform are recommended to be placed in the area marked in blue on the side. It is recommended that the implant sizes to be used should be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental Implants can be used for all bone types. The implants placed in these areas are as follows. For recommended product codes, please refer to MKK.01 and MKK.03 IFU.

Likewise, large diameter implants are also used in areas where the diameter of the tooth and bone is thick, that is, in areas where there are teeth with more than one root. Large-diameter implants are placed in the area of the lower and upper jaw incisors where there is a biting force instead of chewing. In Medisolaris implants; Ø3.7/4.2/4.7/5.2 mm diameter implants for the Merkurhex series with a regular platform and Ø3.7/4.3/5.0 implants for the Venuscon series are recommended to be placed in the area marked in blue on the side. The implant sizes to be used are recommended to be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental Implants can be used for all bone types. The implants placed in these areas are as follows. For recommended product codes, please refer to MKK.01 and MKK.03 IFU.

1.k Bone Type

			
D4	D3	D2	D1
Thin trabecular bone	Thin trabecular bone with thinner porous cortical bone	Dense trabecular bone with porous cortical bone	Dense cortical bone

This surgical guide has been prepared for the various bone densities studied during the use of drills. The dentist is responsible for selecting the protocols to be applied according to the bone quality and anatomical structure during the study.

Bone quality is categorized by bone density. Bone quality is divided into 4 groups as D1, D2, D3 and D4. Bone quality is of great importance in implant dental treatment. It affects many factors such as implant planning, number of implants, surface properties of the implant, implant surgery technique, duration of implant treatment.

D1 bone is the bone quality group in which the bone is the densest and hardest. It is usually seen in the anterior region of the lower jaw, in cases where the bone has melted excessively. It is a group that should be considered in implant treatment. The bone is very hard. Care should be taken not to heat the bone during implant surgery. Otherwise, necrosis of the bone and loss of the implant is inevitable. D1 Bones may have less blood supply, therefore, implant - bone healing will not occur and implants may be lost.

D2 bone is the ideal bone for dental implants. The density and blood supply of the bone and the surface contact between the implant and bone are the best compared to the other groups.

The density of D3 bone is less dense than the density of D2 bone. It is the bone quality that should be approached more cautiously when applying implants.

D4 bone is the bone quality where the density of the bone is close to zero.

In D4 bone quality, a stable fit of the dental implant to the jawbone may not be achieved. In such a case, it is recommended to place the implants in the area to be implanted, after opening the implant socket, compressing implant bone powder into the socket and waiting for 3-6 months.

1.1 Opening of Steril Dental Implant Package

STEP 1: Open the implant box and remove the sterile packaged thermoform.

STEP 2: Open the back of the blister. The product should not be used if the adhesion area is torn, broken or opened. The implant should be opened in a sterile environment where the surgical operation will be performed.

STEP 3: If the implant is to be taken with anguldurva; the key suitable for the anguldurva is attached to the tip of the anguldurva. If the implant is to be taken manually, the manual implant key is inserted into the torque tool.

STEP 4: The cap of the tube is held upwards. After the cap of the tube is turned open, the implant in the stand is removed with the anguldurvaimplant key or manual implant key.

STEP 5: To remove the implant from the stand, the anguldurvaimplant key or manual implant key is inserted into the implant. Make sure that the keys are fully inserted into the implant so that the implant does not fall out.

STEP 6: When the implants are removed from the carrier, they are removed by moving back and forth.

STEP 7: If the anguldurva is used to place the implant, the appropriate torque is adjusted from the motor and placed in the slot opened in the bone. If the implant is placed manually, the implant is removed from the socket and placed in the opened socket after the appropriate torque is adjusted before taking the implant with the torque.

STEP 8: To place a cap on the implant, the cap in the tube is removed in the surgical operation room with the help of a key suitable for the cap.

STEP 9: The cap removed from the tube is placed on the implant without manual contact.

If the products are not to be taken directly with instruments, they should be opened on the sterile cloth.

B. General Characteristics of Drills

All drills are compatible with Venuscon and Merkurhex implants. Drills should be selected according to the diameter and length of the implant to be placed. Among the drills designed for Venuscon and Merkurhex implant system, pilot drills are designed for every size implant with a stopper. With the help of the length lines on the final drills, it is possible to open the appropriate slot for the implant. When the drills are used without a stopper, there are black guide lines on the drills to determine the slot depth. Each of these lines represents a specific length. The lengths they represent are shown in the picture below. The dentist should open the slot using the guide lines according to the length of the implant to be placed. In drills designed for Merkurhex and Venuscon implant system, it is appropriate to use the next drill for hard bones.

All surgical instruments are made of stainless steel.

It is recommended to use each driller maximum 15 times. Before each use, the drills should be cleaned in accordance with the instructions under A.1 General Cleaning Recommendations and sterilized in accordance with the parameters under A.2 Sterilization Process.

Drills should be stored in accordance with A.3 Storage conditions.

B.1 General Cleaning Recommendations

Process Step	Cleaning Process	Applicaition	Equipment and Materials
Step 1	Routine Cleaning	Whenever possible, parts and tools are disassembled before cleaning. Parts are rinsed and wiped with tap water, soaked in warm broad-spectrum enzymatic cleaning solution for at least twenty (20) minutes, and rinsed with tap water.	Wiping cloth, broad spectrum enzymatic cleaning solution and tap water.
Step 2	Disinfection	Using a broad-spectrum enzymatic cleaning solution, the parts are washed, then rinsed and dried.	Wiping cloth, broad spectrum enzymatic cleaning solution and tap water

* As a broad-spectrum enzymatic cleaning solution, it is recommended to use Anios Aniosyme DD1 Concentrated Enzymatic Disinfectant for disinfection.

B.2 Sterilization Process

Surgical Instruments must be sterilized by autoclave method before use. In autoclave sterilization (steam), they should be sterilized at 134°C for 3 minutes.

Method	Steam Sterilization according to ISO 17665
Handover	Pre-Vacum

Temperature	134°C / /270°F
Exposure time	3 Minutes
Pre-Vacuum	3 times < 60 mbar
Minimum Drying time	10 minutes in dry area

The single packaging of Surgical Instruments is not suitable for autoclave sterilization. The products must be unpacked in a controlled environment before sterilization, placed in a package suitable for sterilization and then the sterilization process must be performed.

The surgical set container is suitable for autoclave sterilization. Surgical instruments placed inside this container can be sterilized along with the container.

* Make sure to have the annual maintenance performed on the autoclave to be used for sterilization, and have it calibrated once a year to verify the accuracy of its parameters.

B.3 Storage Conditions

It is recommended that the areas where Surgical Instruments are stored comply with the following temperatures and humidity values.


Temperature	18 C° - 28 C °
Humidity	%40 - %60

B.4 Surgical Instruments

- Punta Drill** : The punta drill is the initial drill and is used to determine where the implant will be placed. There is no coating on the surface. A heat treatment process is applied. The Punto drill is supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.
- Round Drill** : It is the drill used after the punta drill and is used for widening and cleaning the implant placement site after the point. There is no coating on the surface. Heat treatment process is applied. Round drill is supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.
- Merkurhex ve Venuscon Pilot and Final Drills** : The drills are osteotomy drills. Osteotomy drills are used to make a suitable slot in the bone for implant placement. Pilot drills are designed with a stopper to increase safety of use. Pilot drills are 6/8/10/11.5/13/16 mm in length. Pilot drill will be used in soft bones for implants with Ø3.2 diameter. Merkurhex and Venuscon Final Drills are available in five different diameters and 2 lengths. Final drills are designed in 2 types: short and long. Short drills are used where the mouth opening is small and long drills are used where the mouth opening is sufficient. Short drills are used to open slots from 6 mm to 10 mm. Long drills are suitable for slotting up to 16 mm. Drills are designed with two wings. Codes are written on them with laser marking according to the implant diameter to be used. All drills are designed for motorized use. There is no coating on the surface. Heat treatment process is applied. Sharpening process is applied. Drills are supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.
- Countersink Drill** : In cases where the patient's bone is too hard, neck expander drills are used after the final drills to expand the cortical bone and reduce the stress on the neck area of the implant. Countersink drills are available in five different diameters and one size. The countersink drills have coloring and guide lines to guide the diameter and length of the drills. All neck expander drills are designed for motorized use. There is no coating on the surface. Heat treatment process has been applied. Countersink drills are supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.
- Drill Extension** : Used for extending drill lengths. There is no coating on the surface. Heat treatment process is applied. Drill extension is supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.

6. **Rachet Adaptor** : It is an apparatus that allows the use of surgical parts used with a motor. There is no coating on the surface. Heat treatment process is applied. The rachette adaptor is supplied non-sterile and must be sterilized by the end user. Raw material is Stainless Steel.
7. **Transferring and Tightening Part** : Implant Transferring and Tightening Parts are used for placing implants in the implant socket. Short clamping and handling parts are used where the mouth opening is small; long clamping and handling parts are used where the mouth opening is larger. Raw material is Stainless Steel.
8. **Screw Retain Drill**: Screw retain drill are drills used to remove broken implants from the socket in the bone. Raw material is Stainless Steel.
9. **Parallel Pin**: Parallel pins are used to check whether the slots are parallel to each other after placement of multi unit abutments in multiple tooth deficiencies. In multiple implant placement; the parallelism of the implant sockets is important for the upper prosthesis to be made. Straight; 15° and 30° angled types are available. It can be used more than once after cleaning and sterilisation. Raw material is Stainless Steel.
10. **Depth Gauge**: The Depth Gauge is used to check the depth of the socket after the implant sockets are created; before the implant is placed in the socket. After measuring the depth, if the neck is suitable, the implant is inserted; if it is not suitable, the slot is brought to the appropriate size and then the implant is inserted. After cleaning and sterilisation, it can be used more than once. Raw material is Stainless Steel.

The dimensions of the depth gauge are as follows. The unit of measurement is millimeters.

11. **Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool**: It is used so that the Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool, which are used to remove the broken screw in the mouth, can fully centre the screw to be removed. The Broken Screw Removal Centring Part is first placed in the area where the screw to be removed and then the screw is removed from the mouth with the Broken Screw Extraction Drill or / Broken Screw Extractor Hand Tool It is used for removing broken screws in the patient's mouth. After placing the Broken Screw Removal Centring Part, the broken screw is removed with the Broken Screw Extractor Hand Tool. It is released as non-sterile. It can be used maximum 15 times after cleaning and sterilisation. Raw material is Stainless Steel.
12. **Torque Ratchet / Ratchet**: During the placement of implants; during the placement of superstructures with manual parts; It is used for tightening in the mouth according to the recommended torque value. It can be used maximum 15 times after cleaning and sterilisation. Raw material is Stainless Steel.
13. **Hex Drivers** : It is used for carrying implants or abutments into the patient's mouth and tightening them by hand with recommended torques. It is recommended not to use more than 15 times. It can be used up to 15 times after cleaning and sterilisation. Raw material is Stainless Steel.
14. **Handtool** : During the placement of implants in the mouth; It is used with carrying and tightening parts. It can be used maximum 15 times after cleaning and sterilisation. Raw material is Stainless Steel.
15. **Puller Screw** : It is used to remove the abutment stuck in the implant during rehearsal. It can be used after cleaning and sterilisation. Raw material is titanium.
16. **Surgical Sets** : It is a set containing all surgical instruments and instruments used during the preparation of the implant socket, placement of the implant and superstructure. Raw material is Stainless Steel. Surgical sets should only be used in a controlled environment (such as an operating room or an area where appropriate sterilization conditions are maintained).

Products with Measurement Features

Products with Measurement Features	Unit of Measurement under 80/181/EEC	Unit of Measure in the Product Specifications	Tolerance Limit for Products with Measurement Features
Depth Gauge	Milli: m	Mm : Millimeter	Depth Gauge - Lines:

	Meter: m		$\pm 0.03 \text{ mm} / 99,5 \%$
Angled Parallel Pin	Degree: °	Degree: °	Angular Parallel Pin: $\pm 2^\circ$ / 87%

B.5 Preparation of the implant socket

In order to place Medisolaris Dental Implants in the bone, the implant socket should be prepared by following the sequence described in the image below. This protocol applies to soft bones and hard bones. Medisolaris final drills have a single length. In order to open the slot according to the implant size, it is recommended to open the slot by following the depth lines shown by marking after using pilot drills with special stoppers for each size.

To prevent bone heating and burning, it is recommended to use the drills with a physio dispenser that can be supplemented with saline fluid. It should be used after making sure that the drills are fully seated in the angul durwa connected to the physio dispenser. When drilling, drilling should be done using the intermittent drilling technique (pecking method). It is recommended to apply with 700 rpm.

It is recommended to use each driller a maximum of 15 times. Each implant should be considered as one use. Due to the small mouth opening in the molar region, it will be appropriate to use short (DRL xxxS coded) drills and instruments according to the implant to be applied.

A.5.a Recommended Torq Values

Process	Recommended Torq Values
Placement of Implants	35-40 Ncm
Placement of Gingiva Formers	15 Ncm
Placement of Covers	15 Ncm
Tightening of Merkurhex Screw	30 Ncm
Tightening of Venuscon Screw	30 Ncm
Tightening of Multi-Unit Oclusal Screw	20 Ncm
Tightening of Screw Abutments	30 Ncm

The protocol to be followed for the Venuscon implant is as follows.

Venüscon | 3.2mm

3.0/3.2

DRL 2023-16
2.0/2.3mm

Final Drill

DRL 2328
2.3/2.8mm

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CBD 3200

Venüscon | 3.7mm

3.0/3.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

Final Drill

DRL 2834
2.8/3.4mm

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CBD 3700

Venüscon | 4.3mm

3.0/4.3

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

Final Drill

DRL 3438
3.4/3.8mm

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CBD 4200

Venüscon | 5.0mm

3.0/5.0

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

Final Drill

DRL 4146
4.1/4.6mm

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CBD 5000

The protocol to be followed for the Merkurhex implant is as follows.

Merkürhex | 3.2mm

3.0/3.2

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

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CBD 3200

Merkürhex | 3.7mm

3.5/3.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

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CBD 3700

Merkürhex | 4.2mm

3.5/4.2

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

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CBD 4200

Merkürhex | 4.7mm

3.5/4.7

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

DRL 3844
3.8/4.4mm

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CBD 4700

Merkürhex | 5.2mm

3.5/5.2

DRL 2023-16
2.0/2.3mm

DRL 2328
2.3/2.8mm

DRL 2834
2.8/3.4mm

DRL 3438
3.4/3.8mm

DRL 3844
3.8/4.4mm

DRL 4448
4.4/4.8mm

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CBD 5200

B.5 Clinical Benefit: The surgical instruments do not provide a direct therapeutic or diagnostic effect. Their clinical benefit is indirect, as they act as surgical instruments enabling the dentists to prepare bone accurately and safely. The clinical outcome for the patient depends on the dental implant or treatment applied, rather than on the surgical instruments itself.

B.6 Safety and Performance Information: Test reports that meet the safety and performance criteria of Medisolaris products have been applied and the results of the test reports have been found to be appropriate. Negative incidents that may occur in Medisolaris products should be reported to the relevant competent authority.

B.7 Disposal: After the implant removal procedure, the implant component may be infectious, and the instruments may be sharp and piercing, so they are classified as “medical waste.” Implants are classified as “infectious waste”; instruments used in implant placement (e.g., drills, clamping devices) are classified as “sharp waste” due to their cutting and piercing properties.

These containers are manufactured from medium-density polyethylene raw material that is resistant to tearing, puncturing, bursting, and transport, and are leak-proof, double-bottom stitched, and gusset-free, with a double-layer thickness of 100 microns, a minimum lifting capacity of 10 kilograms, and bearing the black "International Biohazard" emblem and the red "CAUTION! MEDICAL WASTE" label in a visible size on both sides. The bags are filled to no more than $\frac{3}{4}$ capacity, tightly sealed, and, where necessary, each bag is placed inside another bag with the same characteristics to ensure complete leak-proofing. These bags cannot be recycled or reused in any way. The contents of medical waste bags must not be compressed in any way; medical waste must not be removed from the bag, emptied, or transferred to another container.

Sharps and piercing medical waste must be collected separately from other medical waste in containers or boxes made of plastic or laminated cardboard with the same properties, bearing the black "International Biohazard" symbol and the words "CAUTION! CUTTING AND PUNCTURING MEDICAL WASTE" in black letters. These collection containers shall be filled to no more than $\frac{3}{4}$ capacity, sealed, and placed in medical waste bags. Once sharps containers are full, they must not be compressed, opened, emptied, or recycled. Implant waste must not be mixed with household waste under any circumstances.

Medical waste collected by hospitals, clinics, dental centers, relevant municipal departments, or companies authorized by municipalities is disposed of by the relevant municipality or authorized medical waste disposal facilities.

Personal protective equipment (gloves, masks, goggles, aprons) must be used before and during disposal. Disposal operations must be carried out in accordance with the Medical Waste Regulation.

Each country must act in accordance with its waste management legislation.

B.8 Meaning of Symbols:

	Batch Code (All Products)		Do not re-use (Dental Implant – Abutment – Gingiva Former – Covers)		Sterilized using irradiation (Dental Implants)
	Catalogue Number (All Products)		Do not re-sterilize (Dental Implant – Abutment – Gingiva Former – Covers)		CE Marking / Notified Body number (All Products)
	Used by date (All Products)		Date of manufacture (All Products)		Non-Steril (All product expect dental implants)
	Do not use if package is damaged and consult instructions for use (All Products)		Caution (All Products)	 (01)01234567891011 (11)200622 (17)240622 (21)88888888	UDI Number (All Products)
	Consult electronic instructions for use (All Products)		Temperature limit (All Products)		Manufacturer (All Products)
	Humidity imitation (All Products)		Keep away from sunlight (All Products)		Double sterile barrier System (Dental Implant)
	Medical Device (All Products)		Unique Device Identifier (All Products)		Keep Dry (All Products)

1-) DENTAL IMPLANTS

Medisolaris (MERKURHEX) Dental Implant and Superstructures produced by MEDİSOLARİS Sağlık Hizmetleri San. Tic. Ltd. Ltd. It is and is surgically placed in the lower or upper jaw bone for the treatment of missing teeth and serves as a replacement for the tooth root by providing a stable basis for restorations. MERKURHEX Dental Implants are bone level dental implants. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material .The target users for Venuscon Dental Implants are dentists.

Superstructures are the device that used to support crowns, bridges or superstructures in the lower and/or upper jaw, used with Dental Implants in completely edentulous or partially edentulous patients.It is recommended that Medisolaris products not be used for purposes other than their stated intended use.

Medisolaris Dental Implant Basic Information : For Merkurhex Series Implant :

Lenght : 6 / 8 / 10 / 11.5 / 13 / 16 mm

Diameter : 3.2/3.7/4.2/4.7/5.2 mm

Lenght Diameter	6 mm	8 mm	10 mm	11.5 mm	13 mm	16 mm
Ø 3.2	Not Available	Available	Available	Available	Available	Available
Ø 3.7	Available	Available	Available	Available	Available	Available
Ø 4.2	Available	Available	Available	Available	Available	Available
Ø 4.7	Available	Available	Available	Available	Available	Available
Ø 5.2	Available	Available	Available	Available	Available	Available

1.a Warning : The information presented here is not sufficient for the immediate application of the Medisolaris Dental Implant System. The user must be legally competent to perform surgical operations and must receive advanced training in the field of oral implantation and detailed knowledge of the system components before use.

1.b Legal Liability Disclosure: The use of this product is the sole responsibility of the end user. No liability is assumed for loss or damage caused by misuse of this product, abnormal use (careless handling, sabotage and ignoring safety information), failure to comply with this IFU. This system can only be used with original Medisolaris Dental parts. The use of third party parts renders Medisolaris Sağlık Hizmetleri A.Ş. products ineffective.

1.c Copyright and Patent: The logos of Medisolaris Dental Implants are protected patents and owned exclusively by Medisolaris Sağlık Hizmetleri A.Ş. any use of these marks is prohibited and is subject to prosecution under Turkish and international law.

1.d Education: Dental implant treatment involves complex surgical and restorative procedures. Such procedures should be performed by dentists who have received the necessary training in this field. Inadequate training can lead to surgical and prosthetic complications, resulting in failure of the implant and restoration.

1.e Patient Population : Those with completed bone development (18 years and above), Not pregnant, not in puberty, not breastfeeding, No hard or soft tissue pathology in the implant area, No active periodontal infection or untreated periodontitis, After periodontal treatment, good oral hygiene, good treatment response, non-smoker, being systemically healthy, No unstable endocrine disease (such as diabetes), No local inflammation or mucosal disease (such as lichen planus), Non-osteoporotic, No immunosuppressive disease, Not requiring cancer treatment (chemotherapy/radiotherapy/anti-resorptive drugs), No acquired or inherited coagulation disorder, Not using/has not used IV-based bisphosphonates. No alcohol or drug addiction, No history of allergy to titanium and stainless steel, Not requiring regular medication due to systemic disease, No mental developmental disorder, Patients without recent myocardial infarction or cerebrovascular accident

1.f Raw Materials and Ingredients: Implants and superstructures are made of ISO 5832-3 Ti6AL4V ELI / GRADE23, surgical instruments are made of 1.4116, hand instruments are made of 1.4057 Stainless Steel. The surface of Medisolaris Dental Implant is roughened by sandblasting and passivation SBM method consisting of hydroxyapatite and tricalcium phosphate to obtain surface microstructure, osseointegration, Ra values 2-2.5Ra.

1.g Indications:

Devices are indicated for total and partially edentulism

1.h Contraindications:

✓ Contraindications ✓ Those whose bone development is not complete (under 18 years of age) ✓ Those whose bone development is not complete (under 18 years of age) ✓ Pregnant, in rapid puberty, or breastfeeding individuals ✓ Presence of hard or soft tissue pathology in the implant area ✓ Individuals with active periodontal infection or untreated periodontitis ✓ Unstable endocrine disease such as uncontrolled diabetes ✓ Presence of local inflammation or mucosal disease such as lichen planus ✓ Individuals diagnosed with osteoporosis	✓ Individuals with immunosuppressive disease ✓ Cancer patients requiring chemotherapy, radiotherapy, or anti-resorptive treatment ✓ Individuals with acquired or hereditary coagulation disorders ✓ Individuals who have used or are using intravenous bisphosphonates ✓ Individuals with alcoholism or substance addiction ✓ Individuals with a history of allergy to titanium or stainless steel ✓ Patients requiring regular medication due to systemic disease ✓ Individuals with mental developmental disorders impairing the ability to maintain oral hygiene ✓ Patients with recent myocardial infarction and cerebrovascular accident
--	---

1.i Residual Risks:


• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability
• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening

1.j Potential Side Effects:

• The patient feels pain/pain	• Occurrence of paresthesia
• Difficulty chewing food	• Bone resorption
• Infection formation,	• Occurrence of bone necrosis
• Spread of infection	• Damage to the adjacent tooth
• Deformation or breakage of the implant or superstructure	• Implant mobilization, loosening, delayed or failure to ossify
• Occurrence of Perforation	• Failure to make a prosthesis
• Tooth loss and need for surgery	• Air embolism
• Pathological condition in the patient	• Numbness, loss of sensation and sensitivity in the patient as a result of nerve damage
• Patient's aesthetic concern	• Loss of stability
• Patient's phonetic anxiety	• Failure to ensure primary stability

• Patient dissatisfaction	• Mandibular fracture
• Excessive bleeding in the patient's implant area	• Hematoma formation
• Obstruction of the patient's airway	• Loss of vitality of adjacent tooth
• Decrease in the patient's quality of life	• Maxillary Sinus Perforation
• Another substance is stuck between the implant and abutment	• Peri implant mukozitis
• Bruising in the implant area	• Postoperative soft tissue complications
• Edema formation in the implant area	• Apart from these, mechanical complications such as nail loosening, cracking, implant cracking, and abutment loosening may also occur.
• Failure (loss) of the implant	• Inability to use or misuse of the product
• Swelling in the implant area	• Allergic reaction
• Exudate in the implant area	• Toxic effect
• Carcinogenic effect	• Wound Opening


1.k Placement Configuration of Implants in the Mouth: Small diameter implants are placed in areas where the alveolar crest is thinner due to the small diameter of the tooth itself and the root. After the implant is placed in the alveolar crest, at least two mm bone distance should remain around the implant in the bucco-lingual and mesio-distal directions. This two mm distance provides the necessary support to the implant during the use of the implant. In Merckurhex series implants; 3.2 mm diameter implants with a narrow platform are recommended to be placed in the area marked in dark colour on the side. The image of the oral structure has been prepared in accordance with the tooth and oral region numbering and nomenclature system defined in ISO 3950. The ISO 3950 tooth and oral region numbering and nomenclature system is shown in the image below. In the image on the right, the areas marked in dark color represent the regions in which the use of implants is recommended according to ISO 3950: teeth numbered 12 and 22 in the maxilla, and teeth numbered 42, 41, 31, and 32 in the mandible. The use of implants is not recommended in individuals with incomplete bone development (under 18 years of age) or in cases where deciduous (primary) teeth are present.

Right					00					Left					Oral cavity									
															01					Maxillary area				
10															20					Quadrant				
03					04										05					Sextant				
18	17	16	15	14	13	12	11	21	22	23	24	25	26	27	28	Permanent teeth								
55					54	53	52	51	61	62	63	64	65					Declduous teeth						
85					84	83	82	81	71	72	73	74	75					Declduous teeth						
48	47	46	45	44	43	42	41	31	32	33	34	35	36	37	38	Permanent teeth								
08					07										06					Sextant				
40															30					Quadrant				
															02					Mandibular area				

It is recommended that the implant lengths to be used should be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME
SLM 3208	Merkurhex Implant Ø3.2mm L:8mm Platform 3.0mm
SLM 3210	Merkurhex Implant Ø3.2 mm L:10 mm Platform 3.0 mm
SLM 3211	Merkurhex Implant Ø3.2 mm L:11.5 mm Platform 3.0 mm
SLM 3213	Merkurhex Implant Ø3.2 mm L:13 mm Platform 3.0 mm
SLM 3216	Merkurhex Implant Ø3.2 mm L:16 mm Platform 3.0 mm

Likewise, large diameter implants are also used in areas where the diameter of the tooth and bone is thick, that is, where there are teeth with more than one root. A large diameter implant is placed in the area where the biting force is located instead of chewing, where the lower and upper jaw incisors are located. In Medisolaris implants, it is recommended to place implants with Ø3.7/4.2/4.7/5.2 mm diameter in the area marked with dark colour on the side for Merckurhex series with regular platform. According to the ISO 3950 tooth and oral region numbering and nomenclature system provided above, the areas marked in dark color in the image on the side represent the regions in which the use of implants is recommended. These regions correspond to teeth numbered 18, 17, 16, 15, 14, 13, 11, 21, 23, 24, 25, 26, 27, and 28 in the maxilla, and teeth numbered 48, 47, 46, 45, 44, 43, 33, 34, 35, 36, 37, and 38 in the mandible. The use of implants is not recommended in individuals with incomplete bone development (under 18 years of age) or in cases where deciduous (primary) teeth are present. The implant sizes to be used are recommended to

be decided by the physician in accordance with the patient's bone structure. Medisolaris Dental implants can be used for all bone types. The implants placed in these areas are as follows.

PRODUCT CODE	PRODUCT NAME	PRODUCT CODE	PRODUCT NAME
SLM 3706	Merkürhex İmplantØ3.7mm L:6mm Platform 3.5mm	SLM 4206	Merkürhex İmplantØ4.2mm L:6mm Platform 3.5mm
SLM 3708	Merkürhex İmplantØ3.7mm L:8mm Platform3.5 mm	SLM 4208	Merkürhex İmplantØ4.2mm L:8mm Platform 3.5mm
SLM 3710	Merkürhex İmplantØ3.7mm L:10mm Platform3.5 mm	SLM 4210	Merkürhex İmplantØ4.2mm L:10mm Platform 3.5mm
SLM 3711	Merkürhex İmplantØ3.7mm L:11.5mm Platform 3.5mm	SLM 4211	Merkürhex İmplantØ4.2mm L:11.5mm Platform 3.5mm
SLM 3713	Merkürhex İmplant Ø3.7 mm L:13 mm Platform 3.5 mm	SLM 4213	Merkürhex İmplant Ø4.2mm L:1 mm Platform 3.5mm
SLM 3716	Merkürhex İmplant Ø3.7 mm L:16 mm Platform 3.5 mm	SLM 4216	Merkürhex İmplant Ø4.2mm L:16mm Platform 3.5mm
SLM 4706	Merkürhex İmplant Ø4.7 mm L:6 mm Platform 3.5 mm	SLM 5206	Merkürhex İmplant Ø5.2mm L:6mm Platform 3.5mm
SLM 4708	Merkürhex İmplant Ø4.7 mm L:8 mm Platform 3.5 mm	SLM 5208	Merkürhex İmplant Ø5.2mm L:8mm Platform 3.5mm
SLM 4710	Merkürhex İmplant Ø4.7 mm L:10 mm Platform 3.5 mm	SLM 5210	Merkürhex İmplant Ø5.2 mm L:10 mm Platform 3.5 mm
SLM 4711	Merkürhex İmplant Ø4.7 mm L:11.5 mm Platform 3.5mm	SLM 5211	Merkürhex İmplant Ø5.2 mm L:11.5 mm Platform 3.5 mm
SLM 4713	Merkürhex İmplant Ø4.7 mm L:13 mm Platform 3.5 mm	SLM 5213	Merkürhex İmplant Ø5.2 mm L:13 mm Platform 3.5 mm
SLM 4716	Merkürhex İmplant Ø4.7 mm L:16 mm Platform 3.5 mm	SLM 5216	Merkürhex İmplant Ø5.2 mm L:16 mm Platform 3.5 mm

1.i Packaging and Sterilization : Medisolaris Dental implants are packaged sterile with Caps. The products are disposable. Gamma radiation is used to sterilize implants after packaging. Thanks to sealed packaging, implants are prevented from being affected by external environments. For proper storage conditions, the temperature should be between 18-28°C and the humidity should be between 40-60%. Sterilization is valid until the expiry time stated on the label on the package. Implants in damaged sterile packages should not be used. Contaminated implants cannot be resterilized. Reuse of implants may result in clinical failure and complications including, but not limited to, infection, failure to osseointegrate implants, bone loss, and damage to adjacent structures. The implant surface should not be touched by hand, as this may deteriorate its surface properties. All other superstructures and hand tools should be sterilized by autoclave. For recommended autoclave sterilization parameters, see the Superstructure section and *MKK.02 Surgical Instruction for Use*. Do not use the product whose packaging is damaged.

Implant tube caps are labeled with diameter markings that are color-coded according to the implant diameter.

Ø 3.2 – Turquase	Ø 3.7 – Green	Ø 4.3 – Purple	Ø 5.0 – Yellow
------------------	---------------	----------------	----------------

1.m Shelf Life: The expiration date is clearly stated on each implant package.

1.n Treatment Planning: Comprehensive patient examination, pre-operative imaging, pre-operative diagnosis and treatment planning should be done before implant treatment. In patients who cannot achieve the required bone density for the implant, it is recommended to apply for bone augmentation procedures with advanced surgical techniques to increase bone density before the application and then perform the application. Inadequate planning and evaluation can lead to implant complications or even implant failure.

1.o Preparation of the Surgical Area: Preparation of the implant area should be done with the surgical drills in the Medisolaris Dental surgical set. The drilling speed should be decided by the physician according to the hardness level of the bone. Drilling should be performed under abundant irrigation. Care should be taken to prevent the bone from getting too hot. Bone temperature should not exceed the upper limit of 47 °C. The use of surgical drills that have lost their sharpness should be avoided. Because these may cause unnecessary trauma to the bone. It is recommended to follow the drilling sequence according to the implant diameter to be used. For detailed information about the surgical operation, see *MKK.02 Surgical Instruction for Use*.

1.p Surgical Placement of Implants: Check the length and diameter of the implant on the product label. Attach the (3) labels included in the packaging to the patient's file. Open the outer cardboard package in a controlled environment and remove the thermoformed packaging. Open the outer cardboard package and remove the thermoformed packaging. Open the tyvek paper on the thermoform and empty the sterile tube into the surgical area. The implant tube consists of a single part: Open the implant tube by turning the top cap. Use the implant carrier in the set to move the implant into the surgically prepared implant slot. The ratchets or physiodispenser included in the set should be used to place the implant. Medisolaris Covers should be closed with the help of a 1.25 mm hex key. Caps and healing caps should be torqued by hand.

1.q Recommended Torq Values

Process	Recommended Torq Values	Process	Recommended Torq Values
Placement of Implants	35-40 Ncm	Screw Venuscon	30 Ncm
Gingiva Former	15 Ncm	Multi-Unit occlusal	20 Ncm
Cover	15 Ncm	Self-Screwing Abutments	30 Ncm

1.r Healing Process of Medisolaris Dental Implants: In general, implants should be left to heal for 2 to 4 months, depending on the patient's bone quality, type and general health status. Immediate loading can be performed with the right case selection and adherence to the immediate loading protocols accepted in the literature. If immediate loading is to be performed, the occlusion environment should be checked. If there are multiple implants, they should be splinted and in cases of complete edentulism, at least 4 implants in the mandible and 6 implants in the maxilla should be connected to each other with a prosthesis. Depending on the bone anatomy, the appropriate number, diameter and length of implants should be placed.

1.s Information to be communicated to the patient after implant placement

It is important that patients follow the instructions given below after dental implant surgery to protect their implants for life. For detailed information, please refer to MKK.04 Patient Education Brochure. It is recommended to be given to patients after the operation.

In dental implant surgical treatment, it is recommended to wait at least 4 months in the upper jaw and at least 3 months in the lower jaw for the implant to fuse with the jaw bone.

After the dental implant is placed, the patient will feel some pain and discomfort. The pain level is usually mild. After the operation, the dentist may prescribe painkillers. With painkillers, the pain can be controlled. These should be taken regularly for 2 days. The first dose of painkillers must be taken before the local anaesthetic wears off. The pain will begin to subside within a few days. If the pain does not go away after a few days, you should consult your doctor.

1. After the dental implant is placed, the patient will feel some pain and discomfort. The pain level is usually mild. After the surgery, the dentist may prescribe painkillers. With painkillers, the pain can be controlled. These should be taken regularly for 2 days. The first dose of painkillers should be taken before the local anesthetic wears off. The pain will start to decrease within a few days. If the pain does not go away after a few days, you should consult your doctor.

2. After the dental implant, the dentist may prescribe antibiotics and antiseptic mouthwash to the patient if deemed appropriate.

3. A sterile sponge is placed on the area where dental implant surgery is performed to control bleeding. The sterile sponge should be bitten by the patient with light force for 20 minutes. At the end of 20 minutes, the sponge should be removed.

4. The patient should not eat or drink anything for two hours after the operation. The patient should not spit or rinse his/her mouth with any liquid for two hours. After two hours, the patient should not eat hard and hot foods on the first day. Hot foods increase both bleeding and the risk of infection.

5. It is normal to have bleeding in the form of leakage during the first 24 hours after the procedure. To control bleeding, the patient bites gauze for 30 minutes. If the bleeding continues for more than 2 days, you should consult our physician. Aspirin and aspirin compounds should not be taken as painkillers because they cause bleeding tendency.

6. Swelling may occur after implant surgery. This is a normal condition. If the pain and swelling do not disappear after 5-6 days, you should consult your physician.

7. The patient should avoid tampering with the area where the dental implant is placed with the tongue, finger or hard instruments.

8. The patient should rest a little after the implant surgery and avoid strenuous exercise for 48 hours.

9. If the patient is using removable dentures, he/she should not use them until the date that the dentist deems appropriate.

10. In dental implant surgery, the wound surfaces are closed with sutures. The sutures must be removed by the dentist after one week at the earliest. The patient is informed by the dentist about the early period in which the treated tooth part can be used.

11. After dental implant surgery, the patient should continue oral care after the first 24 hours. The patient should use antiseptic mouthwash prescribed by the dentist three times a day for the first week after the surgery. Antiseptic mouthwash should not be used after one week. After antiseptic mouthwash, the patient can gargle with warm salt water. Add 1 teaspoon of salt to 1 glass of water and gargle 3-4 times a day to disinfect the area. This will help clean the healing area.

12. After 48 hours, the surgical site should be lightly brushed. Other teeth can be brushed and cleaned normally after the operation. There should be no food residue in the surgical site and the patient should rinse the mouth with warm water after meals.

13. Smoking delays the healing process and prevents dental implants from staying in the mouth for a long time, so the patient should stop smoking cigarettes and similar tobacco before starting dental implant treatment.

14. It is very important to follow all these instructions mentioned above for the success of dental implant surgery. If there are signs of infection after implant surgery or if the patient experiences nausea, fever or vomiting, he/she should consult his/her dentist.

1.t Long-term care of dental implants: In order to be able to use implants for a lifetime without any problems, it is important that they are properly cared for as described above.

1.u The importance of oral hygiene in dental implant patients: Implants are more susceptible to plaque and calculus build-up, inflammation and bone loss than natural teeth. In implant patients, good oral hygiene is very important. The patient should brush teeth correctly and clean the interdental and interimplant areas correctly. Oral hygiene kits are available for the private use of implant patients, which can be used by the patient. If the patient follows the oral hygiene instructions, inflammation of the periimplant tissues will be minimal. Removal of plaque and calculus deposits from the teeth will reduce the risk of periimplantitis, the most common cause of implant failure.

1.v Regular dental visits and examinations for implant patients: In patients treated with dental implant prostheses, check-ups are usually recommended every 3 months for the first 1-2 years and then every 6 months thereafter. At each visit, oral hygiene instructions are reinforced and the patient's ability to maintain good oral hygiene is checked. Patients with implant-retained removable prostheses are given some special instructions. The patient is told how to insert and remove the prosthesis without causing any damage to the implant or prosthesis. Excessive pressure should not be applied to the restoration placed on the implant. Restorations on implants are not as strong as natural teeth. The patient should avoid eating hard foods that can break the restoration. Avoid irritation of the implant area.

1.w Lifespan of dental implants: The lifespan of dental implants depends on many factors, such as patient health and the care provided by the patient. The lifespan of implants depends on how old the patient had the implants, how they live and how they care for them. For the success of the implants, the patient should maintain good oral hygiene and visit the dentist regularly for check-ups. Although a clear lifetime cannot be defined for the device, it is thought to have a lifetime of 10 years under normal conditions. Implants have a shelf life of 5 years and superstructures do not have a shelf life. Implants are an effective solution for missing teeth and can be used for a lifetime with good care. Brushing, flossing and regular visits to the dentist after the implant prosthesis are important to prolong the implant's stay in the mouth. Implants can replace a single tooth, a group of teeth or all upper and lower jaw teeth.

1.x Opening Sterile Packages: To take Medisolaris Dental implants from the product packaging with an angldurva, implant keys suitable for angldurva are used. To take the implants manually, the implant key suitable for manual removal is used. When taking the implants, the size suitable for the area to be placed in the mouth should be selected. For all diameters of the Venuscon implant, the 2.3 mm hexagon crimping spanner included in the surgical set should be used.

In order to attach the closing screws of Medisolaris Dental implants to the implants, the caps that are compatible in diameter that come out of the product are placed with a hex key.

STEP 1: Open the implant box and remove the sterile packaged thermoform.

STEP 2: The back of the thermoform is opened. If the adhesion area is torn, broken or opened, the product should not be used. The implant should be opened in a sterile environment where the surgical operation will be performed.

STEP 3: If the implant is to be taken with an angldurva; the key suitable for the angldurva is attached to the tip of the angldurva. If the implant is to be taken manually, the manual implant key is attached to the torque racket.

STEP 4: The cap of the tube is held upwards. After the cap of the tube is opened by turning the cap, the implant in the sling is removed with an angldurva implant key or manual implant key.

STEP 5: To remove the implant from the sling, the angldurva implant key or manual implant key is inserted into the implant. Make sure that the keys are fully inserted into the implant so that the implant does not fall out.

STEP 6: When the implants are removed from the discs in the sling, they are removed by moving back and forth.

STEP 7: If an angldurva is used to place the implant, the appropriate torque is adjusted from the motor and placed in the slot opened in the bone. If the implant is placed manually, the implant is removed from the slot and placed in the opened slot after the appropriate torque is adjusted before taking the implant to the torque raschette.

STEP 8: In order to place a cover on the implant, the cover in the tube is removed with the help of a key suitable for the cover in the surgical operation room.

STEP 9: The cap removed from the tube is placed on the implant without manual contact.

If the products are not to be taken directly with hand tools, they should be opened on the sterile cloth.

1.y Post-Operative Care: Close follow-up of the patient after surgery is very important. It is recommended to use appropriate oral rinses and to ensure good oral care and not to use tobacco products.

2-) Medisolaris Dental Cover

The closure screw is used to cover and protect the implant after placing the implant in the bone. The closure screw can remain in the mouth for 2-3 months during the surgical healing period.

2.a Raw Material : ISO 5832-23 Ti6AL4V ELI

2.b Indications: The caps come out of the implant package used. It is used in two-stage implant surgery to prevent soft or hard tissue growth on the hollow inner surface of the implant.

2.c Contraindications: Contraindications are mentioned above.

2.d Possible Side Effects: Possible side effects are mentioned above.

2.e Placement: The caps are torqued by hand using a 1.25 mm hex spanner. Recommended torque values can be found in the Recommended Torque Values Table in section **1.p Surgical Placement of Implants**.

2.f Packaging and Sterilisation: The caps are packaged together with the implant in a sterile tube. If the package is intact, sterility is valid until the expiration date indicated on the package. In cases where the cap is packaged individually, the product is single-packed without being sterile and wrapped in aluminum foil; it must be sterilized before use by the autoclave method according to the parameters specified under Section 5, **Sterilization and Packaging**, of the Instructions for Use.

The packaging of the caps is not suitable for autoclave sterilization. Prior to sterilization, the products must be removed from their packaging in a controlled environment, placed into packaging suitable for sterilization, and then sterilized. The calibration of the autoclave must be performed at least once a year, and the results must be within acceptable limits.

3-) Medisolaris Dental Superstructures, Gingiva Formers, Prosthetic Parts and Instruments

Superstructures are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with dental implants in completely edentulous or partially edentulous patients.

3.a Raw Material : ISO 5832-3 / Ti6AL4V ELI (Grade 23) - Stainless Steel 1.4116 and 1.4057

3.b Gingiva Formers: Healing caps are used after the healing of the implant in the jawbone is completed and the overlying soft tissue is removed, until prosthetic treatment for gum shaping. Healing caps remain in the mouth for 15 days or less. They are for Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . single use only.

3.c Temporary Abutment: Temporary abutment is used with Medisolaris implant in completely edentulous or partially edentulous patients. In immediate loading treatments, they are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw for more than 24 hours. Temporary abutment is an endosseous dental implant abutment designed for single use as a temporary prosthesis during the healing process while the permanent prosthesis is being produced. Temporary abutments are compatible with all diameters of Medisolaris implant for Venuscon series. It is fixed to the implant with the help of a screw. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are for single use only.

3.d Straight Abutment: In order to use a cementable abutment, the implant must be placed in the patient without an angle. Cementable Abutments are available for both brands. It is compatible with all diameters of Medisolaris Dental implants for Venuscon series. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is bonded on it with the help of cement. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.e Angled Abutment: Angled abutments facilitate prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid anatomical structure when placing the implant. The angled abutment is compatible with all diameters of the Medisolaris implant for the Venuscon series. It is fixed to the implant with the help of a screw. The angled abutment is used in cemented restorations. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.f Anatomical Straight Abutment: The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This offers a more aesthetic option to the physician. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it emerges through the gum. In a possible gum recession, the underlying titanium is prevented from being seen. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zirconia is glued on it. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is disposable.

3.g Anatomically Angled Abutment : The use of an anatomically angled abutment facilitates prosthesis construction by ensuring that non-aligned sockets are parallel. These abutments help the physician to avoid the anatomical structure when placing the implant. In addition, the use of aesthetic angled abutment shortens the healing time. The anterior and posterior gingival heights of these abutments are different from each other. While the anterior ones are shorter, the posterior ones are longer. This gives the dentist a more aesthetic option. Since the anterior gingival height of the abutment is shorter than the posterior gingival height, the crown on the abutment is deeper; it comes out of the gum. In a possible gum recession, the underlying titanium is prevented from being seen. For each brand (Medisolaris Dental) implant, the same brand Anatomically Angled Abutment should be used. It is fixed to the implant with the help of a screw. A tooth made of materials such as porcelain or zircon is glued on it. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.h Multi Unit Abument: It is the external connection of the implant embedded flat into the bone. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different gum options are available. It is self-screwed and

fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is disposable.

3.i Multi Unit Angled Abutment: It is the external connection of the implant embedded into the bone at an angle. It is used if the physician wants to make the implant screwed over the implant from occlusal. Different angle and gum options are available. It is fixed to the implant with a screw. A tooth is made of materials such as porcelain or zircon. There are two types of multi-unit angled abutment compatible with Medisolaris implants: hex and non-hex. The hex structure is recommended to be used in single tooth restorations, while the non-hex structure is recommended to be used in bridge, bar and hybrid treatments. Hex structure is designed to prevent rotation in single tooth applications. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.j Ball Abutment: It is the external connection of the implant embedded flat into the bone. It is used in cases where the physician wants to apply removable prosthesis to the patient. Different diameters and gum lengths are available. It is self-screwed and fixed to the implant with this screw. A tooth is made of materials such as porcelain or zircon. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.k Tibase Cad-Cam Abutment: It is the platform of the abutment inside the implant. Zirconia or custom (personalised) abutment is produced on it. It is produced for the purpose of using a company-specific platform. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . It is for single use only.

3.l Locater Abutment: It is an abutment system with low vertical height used for implant over implant removable prosthesis in completely edentulous jaws. The locater abutment piece is fixed to the implant with a metal cap and rubber in the removable prosthesis. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.m Universal Straight Abutment: Universal Straight Abutment are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. It is milled and made suitable for the patient. It is fixed to the implant with the help of a screw. When preparing the universal abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant completely or the connection will not be provided. After the universal abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the universal abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. The milled universal abutments should be washed at 55 - 60 °C for 15 minutes. The products should not be sterilised and used without washing. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.n Premill Abutment: It is produced for CAD-CAM devices. The laboratory connects this superstructure to the device and makes a patient-specific abutment. While preparing the premill abutment, the connection area of the abutment should not be prepared. Otherwise, the abutment will not fit the implant properly or the connection will not be provided. After the premill abutment is prepared, vertical grooves (grooves) should be opened on the abutment to prevent the rotation of the crown to be placed on it. Otherwise, the crown may rotate. When preparing the premill abutment, the wall thickness should not fall below 0.20 mm. Otherwise, the mechanical strength of the abutment will decrease. Premilled abutments that have undergone a milling process must be washed at 55–60 °C for 15 minutes. The products must not be sterilized or used without first performing the washing process. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are single-use only.

3.o Ti-base Casting Abutment: For use in personalised cast prosthetic restorations in single or multiple implants. It is fixed to the implant with the help of a screw. Implants for Venuscon series in Medisolaris Dental implants are designed to be compatible with all superstructures. Interface connection type is the same. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.p Screws: The screw allows the superstructures to be mounted to the implant. The screws used with Measuring Parts and Analogues are different. It is packaged and sold together with the related product. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material. Single use only.

3.q Titanium Cap: In completely or partially edentulous patients; These are the parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants. Titanium Caps are designed as Hex and No-Hex. Both caps are compatible with Venuscon Multi Unit Abutments. Ø 5.0 prosthetic platform diameter and 12.0 mm gingival height. No coating process is applied on the surface. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.r Multi-Unit Tibase Abutment : Multi-Unit Tibase abutment; These are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with Venuscon Implants in completely edentulous patients. Multi-Unit Tibase abutment is designed as hex and no-hex. Both abutments are compatible with Venuscon Implants. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . Single use only.

3.s Comfort Cap: The comfort cap prevents the patient from being uncomfortable by covering the multi unit abutment in the patient's mouth until the crown is loaded in Multi unit abutments. Comfort Caps are designed with 5.0 prosthetic platform diameter and 4.8 mm gingival height. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . For single use only.

3.t Metal Ring: Metal rings are parts used to support crowns, bridges or superstructures in the lower and/or upper jaw used with implants in completely edentulous or partially edentulous patients. Metal rings are used in prosthetic restorations together with knob abutment. Metal Rings are designed with 5.0 prosthetic platform diameter and 2.6 mm gingival height. Manufactured from ISO 5832-3 / Ti6AL4V ELI (Grade 23) material . They are for single use only.

3.u Analogs: Analogs are an auxiliary part that enables the patient to restore the chewing function after the Venuscon Dental Implants are surgically placed in the lower or upper jaw bone. Analogs; after the measurements taken with impression pieces are taken on the model; it imitates the implant by mounting on the impression piece. Manufactured from stainless steel 1.4057 material .Single use only.

3.v Impression Part: These are the parts used to determine the dimensions of the crown planned to be placed on the implant and superstructure. There are two separate options for classical and digital methods. Manufactured from Ti6AL4V ELI and stainless steel 1.4057 material . Single use only.

3.w Parallel Pin: In multiple tooth deficiencies; parallel pins are used to check whether the slots are parallel to each other after the placement of multi unit abutments. In multiple implant placement; the parallelism of the implant sockets is important for the upper prosthesis to be made. It can be used more than once after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material

3.x Depth Gauge: The Depth Gauge is used to check the depth of the created slot after the creation of the implant slots; before the implant is placed in the slot. After measuring the depth, if the neck is suitable, the implant is inserted; if not, the slot is brought to the appropriate size and then the implant is inserted. It can be used more than once after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material .

3.y Broken Screw Removal Centering Part / Broken Screw Extractor Hand Tool: It is used to ensure that the Screw Retain Drill and Broken Screw Extractor Hand Tool, which are used to remove the broken screw in the mouth, can centre the screw to be removed. The Broken Screw Removal Centring Part is first placed in the area where the screw to be removed and then the screw is removed from the mouth with the Screw Retain Drill or Broken Screw Extractor Hand Tool. It is used for removing broken screws in the patient's mouth. After placement of the Broken Screw Removal Centering Part, the broken screw is removed with the Broken Screw Extractor Hand Tool. It is released as non-sterile. It can be used maximum 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4057 material

3.z Torque Ratchet / Ratchet : It is used during the placement of implants; during the placement of superstructures with manual parts; for tightening in the mouth according to the recommended torque value. It can be used up to 15 times after cleaning and sterilisation. Manufactured from Ti6AL4V ELI (Grade 23) and Stainless Steel 1.4057 material .

3.aa Hex Drivers: It is used for transporting implants or abutments to the patient's mouth and hand tightening with recommended torques. It is recommended not to use more than 15 times. It can be used maximum 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4116 and stainless steel 1.4057 .

3.bb Hand Tool : During the placement of implants in the mouth; It is used with carrying and clamping parts. It can be used up to 15 times after cleaning and sterilisation. Manufactured from stainless steel 1.4116 and stainless steel 1.4057 .

4. Product Compatibility

For the compatibility of the products with each other and the products recommended for use, please refer to the additional tables. (You can access the tables from the Medisolaris website. www.medisolaris.com) For implant drill compatibility, please refer to MKK.02 Surgical Instruction for Use.

Annex 1 – MKP.02 – L.03 Medisolaris Implant - Superstructures Compatibility List

Annex 2 – MKP.02 – L.04 Medisolaris Implant - Screws Compatibility List

Annex 3 – MKP.02 – L.05 Medisolaris Implant- Surgical Instruments Compatibility List

Annex 4 – MKP.02 – L.13 Medisolaris Implant-Gingival Height/Impression Part/Scan Body/Screwdriver/Transferring and Tightening Part Compliance List

Annex 5 – MKP.02-L.09 Medisolaris Analogue, Impression Part- Screw Compatibility List

Annex 6 – MKP.02-L.10 Medisolaris Multi Unit Abutment - Scan Body -Impression Part - Burn Out Cap Compatibility List

Annex 7 – MKP.02-L.11 Medisolaris Ratchet-Screwdriver/Transferring and Tightening Part Compatibility List

Annex 8 – MKP.02-L.12 Medisolaris Impression Part-Screw Compatibility List

5. Sterilization and Packaging: Superstructure Components, Healing Heads, Prosthetic Parts and Instruments are packaged unsterilised, sealed in aluminium foil. All system parts must be sterilised before use on patients. All parts must be sterilised by autoclave according to the recommended protocols. The recommended steam sterilisation parameters are as in the table below.

Method	Steam Sterilization according to ISO 17665	Exposure time	3 Minutes
Handover	Pre-Vacum	Pre-Vacum	3 times < 60 mbar
Temperature	134°C / /270°F	Minimum Drying Time	10 minutes indoor

The packaging of Superstructure Components, Gingiva Formers, Prosthetic Parts, Auxiliary Parts and Hand Tools are not suitable for autoclave sterilisation. The products must be unpacked in a controlled environment before sterilisation, placed in an autoclave bag and then sterilisation must be performed. The calibration of the autoclave must be performed at least once a year, and the results must be within acceptable limits.

6. Placement: Medisolaris dental abutments are placed on the respective implants with carrying and clamping parts and hex keys. Medisolaris ball abutments are self-screwing and are fixed to the implants in the hex slot above the round part using a ball driver. Recommended torque values for screws and superstructures can be found in the Recommended Torque Values Table in

section **1.p Surgical Placement of Implants.** Cement should not be used in multi-unit abutment application. Ball abutments and Locator abutments should be used for removable prosthesis for patients who are used to removable prosthesis. Angled options of abutments should be used according to the posture of the implant.

7. Clinical Benefit: Dental implants and Surgical Instruments are a safe, functional alternative treatment. For example, it significantly increases the quality of life of the patient, improves the patient's chewing function, greatly reduces the loss of function, and provides relief from symptoms. It is in the form of natural teeth, provides chewing function, provides a better aesthetic appearance, healthy nutrition and self-confidence with healthy teeth.

A strong and stable dental implant restores the missing tooth so that it looks and functions like a natural tooth. Dental implants are advantageous because the stability of removable prostheses can be improved or the need for fixed tooth-supported or removable prostheses can be met. In a study by Perrotti et al., 2 different implant systems were found to be suitable for long-term treatment, with stable marginal bone levels and a 10-year success rate. Dental implants help you maintain the natural shape of your face and smile. Leaving a gap in your mouth after losing one or more teeth leads to further health problems such as jawbone deterioration and loss. When the natural tooth is not supported, the jawbone deteriorates and loses its strength and stability. Dental implants are the only dental restoration option that protects and stimulates natural bone, helps bone to grow and prevents bone loss.

8. Safety and Performance Information: Test reports that meet the safety and performance criteria of Medisolaris products have been applied and the results of the test reports have been found to be appropriate. Adverse events that may occur in Medisolaris products must be reported by the user or patient to the manufacturer and the relevant competent authority.

Performance Information :

Functioning as a tooth root substitute, Designed to eliminate edentulism in the maxilla and mandible, Providing osseointegration in the jawbone by replacing the tooth/teeth

Safety Information :

Please refer to Contraindication, Residual Risks, Warnings, Possible Side Effects in the instructions for use.

9. Warnings:

The causes of tooth loss should be investigated and the continuity of these causes after implantation should be evaluated. Factors that may affect post-implant healing and long-term use include;

► Oral hygiene habits, ► Plaque control efficiency and cognitive skills specific to the individual, ► Periodontal health of existing teeth,

► Soft tissue health, ► Parafunctions.

Patients with the following diseases are recommended to use implants after being checked by physicians (with consultation if necessary);

1. In patients with heart disease (those using a pacemaker), implant treatment is recommended only after a cardiology consultation. (Subject to the physician's evaluation.)
2. In patients with hypertension, blood pressure should be measured and evaluated before implant treatment. The patient should be informed about the regular use of antihypertensive medications. (Subject to the physician's additional evaluation.)
3. Patients who have not been diagnosed with osteoporosis but are at risk: implant treatment should be performed after an osteoporosis consultation. (Subject to the physician's additional evaluation.)
4. In patients who have used oral medications for more than 3 years, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
5. In patients with immune system disorders or immunosuppressed patients, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
6. In patients who have previously received radiotherapy or chemotherapy, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
7. In patients with psychological and mental disorders, implant treatment should be performed after consultation. (Subject to the physician's additional evaluation.)
8. Patients who smoke, consume alcohol, or use medications; and patients with bruxism (teeth clenching or grinding).

The implant site must be positioned correctly.

The distance from adjacent natural teeth must be at least 1.5 mm, and the distance from another dental implant must be at

least 3 mm.

The amount of bone on the vestibular, lingual, or palatal side must be at least 1 mm.

During preparation of the implant bed, the bone temperature must not exceed 47 °C.

Surgical instruments must be sterilized before use (according to the manufacturer's recommendations).

Surgical instruments must not be used beyond their maximum number of uses.

The user must be legally qualified to perform surgical procedures.

Before use, detailed information about oral implantology and system components must be obtained through advanced training.

The patient must be provided with a Patient Information Brochure.

For pain management, the prescription written by the dentist should be used regularly for the duration specified by the physician.

If the patient's pain does not subside within a few days, they must consult the physician.

The sterile gauze placed at the surgical site for bleeding control after dental implant surgery should be bitten with light pressure for 20 minutes, unless otherwise instructed by the physician, and removed after 20 minutes.

For two hours after surgery, the patient should not eat or drink anything. During these two hours (unless otherwise instructed by the physician), the patient should not spit or rinse the mouth with any liquid. For the rest of the first day (after the initial two hours), hard and hot foods should be avoided.

The implant site must not be touched with the tongue, fingers, or any hard object.

After the operation, heavy sports should be avoided and the patient should rest.

An antiseptic mouthwash should be used for one week after implant treatment.

Tobacco products must not be used after implant treatment.

Implants in damaged sterile packages must not be used.

Contaminated implants must not be re-sterilized.

Reusing implants may lead to clinical failure and complications including, but not limited to, infection, lack of implant integration, bone loss, and damage to adjacent structures.

A comprehensive patient examination, preoperative imaging, preoperative diagnosis, and treatment planning must be performed before implant treatment. Inadequate planning and evaluation may lead to implant complications or even implant failure.

MR-Conditional

Non-clinical testing has demonstrated that implants belonging to the Medisolaris Venuscon and Merkurhex product families, when used in combination with straight and angled abutments, can be conditionally scanned in whole-body MRI systems with a 3-Tesla static magnetic field. A patient with these devices may undergo an MRI scan under the following conditions:

- **Static Magnetic Field:** 3-Tesla
- **Maximum SAR Value:** Whole-body averaged Specific Absorption Rate (SAR) reported by the scanner of 2 W/kg
- **Magnetically Induced Displacement Angle:** 2°
- **Scan Duration:** 15 minutes of continuous scanning
- **RF-Induced Heating:** Under the scanning conditions specified above, the maximum temperature increase generated by the devices in the worst-case scenario is 2.6°C for the Venuscon family and 2.4°C for the Merkurhex family
- **Image Artifact:** The quality of the MR image may be affected by the presence of the implant. In testing, the maximum artifact extent was measured as 11 mm for Spin Echo sequences and 19 mm for Gradient Echo sequences. Therefore, if the imaging region is close to the implant location, optimization of MRI parameters is recommended
- **Mechanical Safety:** It has been verified that the magnetically induced displacement force and torque are lower than the gravitational forces encountered during normal daily activities and do not pose a clinical risk

Caution: This safety information applies only to the tested Medisolaris implant and abutment combinations; MR safety cannot be guaranteed if the device is used in combination with components from other manufacturers.

10. Information to be given to the patient : MKK.04 Patient Education Brochure must be given to patients by physicians.

11. Disposal: After the implant removal procedure, the implant component may be infectious, and the instruments may be sharp and piercing, so they are classified as "medical waste." Implants are classified as "infectious waste"; instruments used in implant placement (e.g., drills, clamping devices) are classified as "sharp waste" due to their cutting and piercing properties.

These containers are manufactured from medium-density polyethylene raw material that is resistant to tearing, puncturing, bursting, and transport, and are leak-proof, double-bottom stitched, and gusset-free, with a double-layer thickness of 100 microns, a minimum lifting capacity of 10 kilograms, and bearing the black "International Biohazard" emblem and the red "CAUTION! MEDICAL WASTE" label in a visible size on both sides. The bags are filled to no more than ¾ capacity, tightly sealed, and, where necessary, each bag is placed inside another bag with the same characteristics to ensure complete leak-

proofing. These bags cannot be recycled or reused in any way. The contents of medical waste bags must not be compressed in any way; medical waste must not be removed from the bag, emptied, or transferred to another container.

Sharps and piercing medical waste must be collected separately from other medical waste in containers or boxes made of plastic or laminated cardboard with the same properties, bearing the black "International Biohazard" symbol and the words "CAUTION! CUTTING AND PUNCTURING MEDICAL WASTE" in black letters. These collection containers shall be filled to no more than ¾ capacity, sealed, and placed in medical waste bags. Once sharps containers are full, they must not be compressed, opened, emptied, or recycled. Implant waste must not be mixed with household waste under any circumstances.

Medical waste collected by hospitals, clinics, dental centers, relevant municipal departments, or companies authorized by municipalities is disposed of by the relevant municipality or authorized medical waste disposal facilities.


Personal protective equipment (gloves, masks, goggles, aprons) must be used before and during disposal. Disposal operations must be carried out in accordance with the Medical Waste Regulation.

Each country must act in accordance with its waste management legislation.

12. SCCP

Following the certification of Medisolaris products under MDR 2017/745, the SCCP link will be added.

13. MEANING OF SYMBOLS

	Batch Code (All Products)		Do not re-use (Dental Implant – Abutment – Gingiva Former – Covers)		Sterilized using irradiation (Dental Implants)
	Catalogue Number (All Products)		Do not re-sterilize (Dental Implant – Abutment – Gingiva Former – Covers)		CE Marking / Notified Body number (All Products)
	Used by date (All Products)		Date of manufacture (All Products)		Non-Steril (All product expect dental implants)
	Do not use if package is damaged and consult instructions for use (All Products)		Caution (All Products)		UDI Number (All Products)
	Consult electronic instructions for use (All Products)		Temperature limit (All Products)		Manufacturer (All Products)
	Humidity imitation (All Products)		Keep away from sunlight (All Products)		Double sterile barrier System (Dental Implant)
	Medical Device (All Products)		Unique Device Identifier (All Products)		Keep Dry (All Products)
	MRI Conditional (Dental Implant – Abutment – Screw)				

MEDİSOLARİS SAĞLIK HİZMETLERİ SANAYİ ve TİCARET LİMİTED ŞİRKETİ - Fatih Mahallesi 1191. Sokak No: 28/1 B
Blok Apt. Gaziemir / İZMİR - +90 232 203 6655 / www.medisolaris.com