

İMPLANT-CERRAHİ EL ALETLERİ UYUM LİSTESİ

Hazırlayan : Hande AŞAR
Yönetim Temsilcisi

Onaylayan : Basri EROL
Genel Müdür

Hazırlanma Tarihi : 12.09.2020
Rev No : 001
Rev Tarihi : 12.10.2021

DRİLL	İMPLANT																							
	SLV 3208	SLV 3210	SLV 3211	SLV 3213	SLV 3216	SLV 3706	SLV 3708	SLV 3710	SLV 3711	SLV 3713	SLV 3716	SLV 4306	SLV 4308	SLV 4310	SLV 4311	SLV 4313	SLV 4316	SLV 5006	SLV 5008	SLV 5010	SLV 5011	SLV 5013	SLV 5016	
DRL 2023-6	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Red
DRL 2023-8	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red
DRL 2023-10	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red
DRL 2023-11	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Red	Green	Red	Red	Red
DRL 2023-13	Red	Red	Red	Green	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Red	Red	Green	Red
DRL 2023-16	Red	Red	Red	Red	Red	Red	Red	Red	Red	Red	Green	Red	Red	Red	Red	Green	Red	Red	Red	Red	Red	Red	Red	Green
DRL 2328S	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Red
DRL 2328	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green
DRL 2834S	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Red
DRL 2834	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green
DRL 3438S	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Red
DRL 3438	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green
DRL 3844S	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Red
DRL 3844	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green	Green
DRL 4146S	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Green	Green	Green	Red	Red	Red	Red

